

Corning® FluoroBlok™ Inserts

Frequently Asked Questions

CORNING

The Corning FluoroBlok membrane is black and has improved spectral characteristics. General information in this FAQ applies to both black and former purple versions, but specific wavelength ranges apply to the black version.

For details, see PET Membrane for Corning FluoroBlok 3.0 µm and 8.0 µm Pore Size Cell Culture Inserts (CLS-DL-CC-042).

Corning FluoroBlok inserts have a black dyed polyethylene terephthalate (PET) microporous membrane that blocks light transmission at visible wavelengths between 400 to 700 nm. Migration assays are performed in the traditional manner using this insert system, with cells that are labeled via a fluorescent dye. This allows researchers to specifically detect and quantify fluorescently labeled cells that have migrated through the insert using a fluorescence plate reader or an inverted fluorescence microscope.

Corning FluoroBlok inserts can be used to study a wide range of cell types and activities such as:

- ▶ Inflammation with neutrophils,¹⁻⁵ transepithelial⁶ and transendothelial⁷ migration; and analysis of blood-brain barrier,^{8,9} dendritic cells,¹⁰ and Macrophages¹¹
- ▶ Pathways for stem cell differentiation^{12,13}
- ▶ Screening for population-specific neuronal motogens¹⁴
- ▶ Migration of normal, transformed and transfected cells^{15,16}
- ▶ Chemoinvasion assays, drug screening^{17, 18}

1. Can I do tumor cell invasion with Corning FluoroBlok inserts?

Yes. You can coat the FluoroBlok inserts with Corning Basement Membrane matrix, just as you would a clear insert.

Visit www.corning.com/lifesciences to access a protocol for Cell Invasion Assay (CLS-DL-CC-031) as a guide for self-coating the membrane.

2. Can I study angiogenesis cell migration and invasion with the Corning FluoroBlok inserts?

Yes. The following product can be used to study angiogenesis cell migration using FluoroBlok inserts.

Corning BioCoat Angiogenesis System: Endothelial Cell Migration 3.0 µm





Cat. No.	Description	Qty/Cs
354144	Five insert plates with 24-well plates and lids	5

3. Does the Corning FluoroBlok membrane autofluoresce?

The FluoroBlok membrane exhibits negligible autofluorescence across the visible spectrum (400 to 700 nm) as demonstrated by top-reading fluorescence data. However, there is a low level of fluorescence background in bottom reading mode due to autofluorescence of and/or reflection from the polystyrene well bottom of the base plate. Use of excessively high gain settings or failure to run the appropriate controls can often give the false impression that the FluoroBlok membrane blocks light inefficiently or has high inherent autofluorescence.

Tumor Invasion Model Assay

HT-1080 and 3T3 cells post-labeled with Calcein AM after migration through uncoated 8.0 μm inserts (Figures A and C) and invasion through Corning Matrigel[®] matrix coated 8.0 μm inserts (Figures B and D). HT-1080 cells are capable of migration (A) and invasion (B). Lacking matrix metalloproteinases (MMPs), 3T3 cells can migrate (C) but did not invade well (D) through the basement membrane.



4. What is the advantage of using a Corning[®] FluoroBlok™ insert instead of a typical clear cell insert?

FluoroBlok inserts allow researchers to detect fluorescently labeled cells, or cells transfected to express fluorescent proteins, passing through the membrane in a homogeneous format (i.e., no further cell separation, washing or harvesting is necessary to detect cells specifically passing through the membrane). Simply add cells, stain (if necessary), and read cells on a bottom-reading plate reader.

5. Why would I want to use a homogeneous assay with Corning FluoroBlok inserts?

Many cell migration and invasion assays require destructive, time-consuming manual processing of inserts for cell quantitation. Examples include using clear inserts for cell migration studies (cells must be removed from the lower well and labeled in some way for analysis) and invasion studies (cells are removed from the top side of the membrane with a cotton swab, and cells on the underside are then counted). With FluoroBlok inserts, fluorescently labeled cells can be quantified directly and in real time with a fluorescence plate reader by reading the desired well for fluorescence output. For invasion and migration assays, no tedious manual sample processing is necessary, so assays may be automated for high throughput screening.

6. How do Corning FluoroBlok inserts work?

The membrane within the insert blocks light transmission from 400 to 700 nm. By labeling cells with a fluorophore that has either excitation or emission wavelength within this range, it is possible to specifically quantify migration or invasion below the membrane. This can be done by exciting the well from below while simultaneously measuring fluorescence emission from below the membrane. The membrane prevents cells above the membrane from becoming excited or emitting and influencing the signal measured from below.

7. I notice that when I look through the Corning FluoroBlok membrane, it does not appear to be opaque. How can it block light if I can see through it?

You are seeing light passing through the pores in the membrane. This is normal.

8. What types of fluorescent dyes can I use to label cells and detect them with the Corning FluoroBlok insert?

Any dye that has an emission wavelength between 400 and 700 nm can be used with this system with high confidence. Outside of this spectral range, there will be more of a contribution from cells remaining in the apical chamber. Corning Calcein AM fluorescent dye (Cat. Nos. 354216 and 354217) and Corning DiIC₁₂(3) fluorescent dye (Cat. No. 354218), are available for labeling cells. For information about these and other compatible dyes, visit www.corning.com/lifesciences.

9. Are the Corning FluoroBlok inserts compatible with any plate?

No. It is critical to use the proper Falcon[®] receiver plate with the proper Corning FluoroBlok insert system. Individual inserts must be used with Falcon cell culture insert companion plates (24-well, Cat. No. 353504) to properly position the inserts in the wells. Corning FluoroBlok 24 Multiwell inserts use the supplied Falcon 24-well plates. These plates accurately fit the Corning FluoroBlok 24 Multiwell insert. The Corning FluoroBlok 96 Multiwell inserts use the supplied square well receiver plates. These plates must be used for running the assay, labeling, and reading samples to achieve reliable assay results.

10. Can Corning FluoroBlok inserts be coated?

If desired, the inserts can be coated using Corning ECM proteins, such as Corning Matrigel[®] matrix for cell invasion, or Fibronectin for migration. Uncoated Corning FluoroBlok inserts are available in a variety of formats and pore sizes.

Products may not be available in all markets.

Corning FluoroBlok Cell Culture Inserts

Cat. No.	Description	Qty/CS
351151	3.0 μm inserts for 24-well plates	48
351152	8.0 μm inserts for 24-well plates	48

Corning FluoroBlok 24-Multiwell Insert Systems, 3.0 μm

Cat. No.	Description	Qty/CS
351156	Five insert plates with 24-well plates and lids	5
Corning FluoroBlok 24-Multiwell Insert Systems, 8.0 μm		
351157	One insert plate with 24-well plate and lid	1
351158	Five insert plates with 24-well plates and lids	5
Corning FluoroBlok 96-Multiwell Insert Systems, 3.0 μm		
351161	One insert plate with 96-well plate and lid	1
351162	Five insert plates with 96-well plates and lids	5
Corning FluoroBlok 96-Multiwell Insert Systems, 8.0 μm		
351163	One insert plate with 96-well plate and lid	1
351164	Five insert plates with 96-well plates and lids	5

For assistance selecting a pre-coated insert or ECM for your application, or coating your Corning FluoroBlok inserts, contact ScientificSupport@corning.com.

11. How can I be sure that the dye won't leach out of the membrane and contaminate my sample?

Corning® FluoroBlok™ (black) membranes and insert systems were tested using commonly used biological solvents including: saline, 10% DMSO in culture medium, 4% paraformaldehyde, and 100% methanol. The effects on total transmission were negligible, and no observable dye was leached from the membranes.

12. When I set up my plate reader, can I use standard 24-well and 96-well templates?

No. It is critical that the detector is properly positioned under the Corning FluoroBlok inserts. Visit www.corning.com/lifesciences for proper set up of your plate reader. If your plate reader is not listed, contact ScientificSupport@corning.com.

13. Can I use Corning FluoroBlok inserts for the migration/invasion of suspension cells such as lymphocytes?

Yes. The inserts work well for non-adherent quickly migrating cells. If you pre-label your cells, you can collect kinetic data. Visit www.corning.com/lifesciences and the references listed.¹⁻¹¹

14. Can cells be removed from the inserts after migration/invasion?

Yes. If the cells are adherent, use of trypsin, Corning Dispase (Cat. No. 354235), Corning Cell Recovery Solution (Cat. No. 354253), Accutase® (Cat. No. 25-058-CI) or other enzymatic methods may prove successful. Refer to Corning Matrigel® Matrix Frequently Asked Questions (CLS-DL-CC-026) at www.corning.com/lifesciences to see which cell dissociation reagents to use based on downstream applications.

References

- Nick JA, et al. Recombinant human activated protein C reduces human endotoxin-induced pulmonary inflammation via inhibition of neutrophil chemotaxis. *Blood* 104(13):3878 (2004).
- Kuijpers TW, et al. Apoptotic neutrophils in the circulation of patients with glycogen storage disease type 1b (GSD1b). *Blood* 101(12):5021 (2003).
- Kuijpers TW, et al. Neutrophils in Barth syndrome (BTHS) avidly bind annexin-V in the absence of apoptosis. *Blood* 103(10):3915 (2004).
- Slofstra SH, et al. Inhalation of activated protein C inhibits endotoxin-induced pulmonary inflammation in mice independent of neutrophil recruitment. *Br J Pharmacol.* 149(6):740 (2006).
- Stegenga ME, et al. Effect of acute hyperglycaemia and/or hyperinsulinaemia on proinflammatory gene expression, cytokine production and neutrophil function in humans. *Diabet Med.* 25(2):157 (2008).
- Lin A, et al. Streptolysin S Inhibits Neutrophil Recruitment during the Early Stages of Streptococcus pyogenes. *Infect Immun.* 77(11):5190 (2009).
- Johnson LA, et al. An inflammation-induced mechanism for leukocyte transmigration across lymphatic vessel endothelium. *JEM* 203(12):2763 (2006).
- Chaudhuri A, STAT1 signaling modulates HIV-1-induced inflammatory responses and leukocyte transmigration across the blood-brain barrier. *Blood* 111(4):2062 (2008).
- Chaudhuri A, 1 HIV-1 activates proinflammatory and interferon-inducible genes in human brain microvascular

- endothelial cells: putative mechanisms of blood-brain barrier dysfunction. *J Cereb Blood Flow Metab.* 28:697 (2008).
10. Lapteva N, et al. Attraction and activation of dendritic cells at the site of tumor elicits potent antitumor Immunity. *Molec Ther.* 17(9):1626 (2009).
 11. Link TM, et al. TRPV2 has a pivotal role in macrophage particle binding and phagocytosis. *Nat Immun.* 11:232 (2010).
 12. Brown MD, et al. Influence of omega-6 PUFA arachidonic acid and bone marrow adipocytes on metastatic spread from prostate cancer. *Br. J. Cancer.* 102:403 (2010).
 13. De Becker A, et al. Migration of culture-expanded human mesenchymal stem cells through bone marrow endothelium is regulated by matrix metalloproteinase-2 and tissue inhibitor of metalloproteinase-3. *Haematologica* 92(4):440 (2007).
 14. Hassoun AT, et al. A rapid screening method for population-specific neuronal motogens, substrates and associated signaling pathways. *J. Neurosci. Methods.* 166:178 (2007).
 15. Bernhagen J, et al. MIF is a noncognate ligand of CXC chemokine receptors in inflammatory and atherogenic cell recruitment. *Nat Med.* 13:587 (2007).
 16. Vitari AC, et al. COP1 is a tumour suppressor that causes degradation of ETS transcription factors. *Nature* Published online 15 May 2011.
 17. Albini A, Noonan DM. The “chemoinvasion” assay, 25 years and still going strong: the use of reconstituted basement membranes to study cell invasion and angiogenesis. *Curr Opin Cell Biol.* 22(5):677-89 (2010).
 18. Payen VL, et al. Monocarboxylate Transporter MCT1 promotes tumor metastasis independently of its activity as a lactate transporter. *Cancer Res.* 77(20):5591-5601 (2017).

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences
www.corning.com/lifesciences

NORTH AMERICA
t 800.492.1110
t 978.442.2200

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India
t 91 124 4604000
Japan
t 81 3-3586 1996
Korea
t 82 2-796-9500
Singapore
t 65 6572-9740
Taiwan
t 886 2-2716-0338

EUROPE
CSEurope@corning.com
France
t 0800 916 882
Germany
t 0800 101 1153
The Netherlands
t 020 655 79 28
United Kingdom
t 0800 376 8660

All Other European Countries
t +31 (0) 206 59 60 51

LATIN AMERICA
grupoLA@corning.com
Brazil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400