

# Corning® BioCoat™ High Content Imaging Glass Bottom Microplates for High Throughput Data Capture and Analysis

CORNING

## Technical Brief

Hilary Sherman  
Corning Incorporated  
Life Sciences  
Kennebunkport, ME 03820

Corning® BioCoat™ high content imaging glass bottom microplates are ideal for performing cell-based imaging assays where cell attachment or retention may be an issue with a standard glass microplate. The BioCoat high content imaging microplates are available pre-coated with rat tail collagen type I, human fibronectin, or poly-D-lysine.

### Benefits

- ▶ Uniform and consistent coating for improved cell attachment and even distribution, or spreading of cells
- ▶ High optical clarity and scratch-resistant glass
- ▶ Glass bottom thickness of 200  $\mu\text{m}$  is well-suited for imaging microscopy
- ▶ Well bottom flatness <50  $\mu\text{m}$  to ensure planarity for imaging devices
- ▶ Low background fluorescence and minimal crosstalk to provide outstanding optical quality for cell-based assays

### Application

Glass bottom microplates have historically been viewed as the gold standard for high content cell-based assays because of the high optical clarity, flatness, and scratch resistant properties intrinsic to glass. Conversely, the natural hydrophobicity of glass may cause adherence issues with some cell types. In those cases, cells may benefit greatly from a protein coating to aid in cell attachment, spreading, and/or functionality. Studies described here demonstrate that Corning BioCoat high content imaging glass bottom microplates provide improved cell attachment, spreading, and retention.

HepG2 cells were seeded on uncoated and collagen-coated high content imaging glass bottom microplates. Cells cultured on the uncoated glass microplates were clumpy and displayed an uneven distribution compared to cells cultured on the collagen-coated microplates. (Figures 1 and 2). As a result, cells cultured on the uncoated plates exhibited statistically lower cell counts when analyzed using a high content imager (data not shown).

In a separate experiment, HEK-293 cells were seeded on fibronectin and collagen-coated glass microplates. After multiple washes, a greater number of cells were retained on fibronectin and collagen-coated glass microplates when compared to the uncoated microplates (Figure 3).

For more details on this experiment and other assays performed using Corning BioCoat high content imaging microplates, please refer to Corning document CLS-AN-244.




Figure 1. HepG2 cells on uncoated glass microplates were clumpy and displayed an uneven distribution.




Figure 2. HepG2 cells on coated glass microplates displayed an even distribution.





Figure 3. HEK-293 cells exhibit improved retention on Corning BioCoat Fibronectin and Collagen I high content imaging glass bottom microplates compared to uncoated glass microplates. Data shown with standard errors. One-way ANOVA with Newman-Keuls post-test \*\*\* $p < 0.001$ .  $n = 768$  from 2 independent studies. 16 fields per well were analyzed.

## Ordering Information

### Corning® BioCoat™ High Content Imaging Glass Bottom Microplates

| Cat. No. | Description  | Sterile | Qty/Pk | Qty/Cs |
|----------|--|---------|--------|--------|
| 4582     | 96-well half area microplate, glass bottom, collagen-coated, with lid      | No      | 1      | 10     |
| 4583     | 384-well microplate, glass bottom, collagen-coated, with lid               | No      | 1      | 10     |
| 4584     | 96-well half area microplate, glass bottom, fibronectin-coated, with lid   | No      | 1      | 10     |
| 4585     | 384-well microplate, glass bottom, fibronectin-coated, with lid            | No      | 1      | 10     |
| 4586     | 96-well half area microplate, glass bottom, poly-D-lysine-coated, with lid | No      | 1      | 10     |
| 4587     | 384-well microplate, glass bottom, poly-D-lysine-coated, with lid          | No      | 1      | 10     |

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.



At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please call 800.492.1110 or visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences). Customers outside the United States, call +1.978.442.2200 or contact your local Corning sales office listed below.

#### Corning Incorporated Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

#### Worldwide Support Offices

##### ASIA/PACIFIC

**Australia/New Zealand**  
t 0402-794-347

**China**  
t 86 21 2215 2888  
f 86 21 6215 2988

**India**  
t 91 124 4604000  
f 91 124 4604099

##### Japan

t 81 3-3586 1996  
f 81 3-3586 1291

##### Korea

t 82 2-796-9500  
f 82 2-796-9300

##### Singapore

t 65 6733-6511  
f 65 6861-2913

##### Taiwan

t 886 2-2716-0338  
f 886 2-2516-7500

##### EUROPE

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

##### All Other European Countries

t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

##### LATIN AMERICA

**Brasil**  
t (55-11) 3089-7419  
f (55-11) 3167-0700

**Mexico**  
t (52-81) 8158-8400  
f (52-81) 8313-8589

**CORNING** | **FALCON** | **AXYGEN** | **GOSELIN** | **PYREX**

For a listing of trademarks, visit us at [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks). All other trademarks in this document are the property of their respective owners.