

Corning® Elplasia® Plates

Assays Requiring Multiple Spheroids per Well

Application Note

Hilary Sherman and John Shyu
Corning Incorporated, Life Sciences
Kennebunk, ME USA

CORNING

Introduction

It has been well-established that 3D cancer models are more reflective of *in vivo* tumor microenvironment when compared to traditional 2D models¹. As a result, there is a growing trend for more representative 3D cancer cell models. Traditional scaffold-free 3D models often force researchers to choose between a) multiple spheroids per well of a non-uniform shape and size, or b) a single uniform spheroid per well. However, multiple uniform spheroids in a single well are often desired to increase assay signal or increase the number of data points under a single condition. Here we demonstrate how Corning Elplasia 96-well microplates can reproducibly form multiple uniform ovarian cancer spheroids that can be used for luminescent, fluorescent, and high content imaging assays.

Materials and Methods

Spheroid Formation

Corning Elplasia plates (Corning 4442) were pre-wet prior to seeding cells by adding 50 μ L of DMEM (Corning 10-013-CM) with 10% fetal bovine serum (Corning 35-010-CV) per well and centrifuging at 500 x g for 1 min. Once all trapped air was removed, HEY-T30 (ATCC® CRL-3252™) or SKOV3 (ATCC HTB-77) cells were seeded into plates at approximately 100 cells per microcavity in a volume of 150 μ L per well. Cells were incubated for 24 hours.

Spheroid Consistency

After an overnight incubation of HEY-T30 spheroid formation, 20 μ L of NucRed® Live 647 ReadyProbes® Reagent (Thermo Fisher R37106) was added to each well and incubated for 20 to 30 minutes. Once cells were completely stained, spheroids were imaged with the Thermo Scientific CellInsight™ CX7 confocal imager to assess spheroid circularity and area.

Homogenous Cell-based Assay

After an overnight incubation for HEY-T30 spheroid formation, 10 μ L of Triton™ X solution to achieve a final concentration of 0.5% or buffer was added into each well. Sixteen to 24 hours later, equal volume of CellTiter-GLO® 3D Cell Viability Assay (Promega G9683) was added to each well and allowed to incubate at room temperature for 30 minutes with a 5-minute shake step. Plates were then read for luminescent signal using a PerkinElmer Envision.

High Content Imaging Assay

After an overnight incubation, HEY-T30 and SKOV3 spheroids were exposed to cisplatin (Thermo Fisher 1134357) by adding 20 μ L at various dilutions. The plates were incubated overnight. In order to assess cytotoxicity, 10 μ L per well of 340 μ g/mL Hoechst 34580 (Thermo Fisher H21486) and 136 μ g/mL propidium iodide (AnaSpec 83215) diluted in phosphate buffered saline (PBS; Corning 21-030-CM) were added to each well for 1 hour. Plates were then imaged on the Thermo Scientific CellInsight™ CX7 confocal imager.

Results and Discussion

Assay Robustness

The unique design of Corning Elplasia plates promote the generation of multiple, uniform spheroids in each well (Figure 1). To assess the ability of Corning Elplasia plates to generate uniform spheroids, HEY-T30 cells were seeded into plates overnight. Once spheroids were generated, they were stained with a fluorescent dye so that circularity and consistency of the spheroids could be examined (Figure 2). The data shows that the spheroids generated were consistent in shape and size throughout the plate. In addition to consistency of spheroid generation, it is also important to be able to execute homogenous cell-based assays in Corning Elplasia plates. Figures 3 and 4 demonstrate high signal-to-background ratios (S/B) and repeatable Z' from a homogenous luminescent assay run directly in the plate.




Figure 1. Uniform, single spheroids formed in each microcavity. Representative image of one field from each well of a 96-well Corning Elplasia plate with one field digitally zoomed in. Images were taken with the Thermo Scientific CellInsight CX7 confocal imager using a 4X objective.




Figure 2. Uniform spheroids generated in each microcavity. Representative image of NucRed stained HEY-T30 spheroids (A) used to assess spheroid circularity (B) and average spheroid area per well (C) from several studies. Images were taken with the Thermo Scientific CellInsight CX7 confocal imager using a 4X objective. Data shown with standard deviation from 3 independent studies. N=288 wells. Scale bar is 100 μm.




Figure 3. Large luminescence assay window. Signal-to-background ratios from 3 independent tests using HEY-T30 spheroids exposed to 0.5% Triton X or buffer.




Figure 4. Excellent Z'. Z' values above 0.7 were achieved with 3 independent studies.

Proof of Concept Assay

To demonstrate Corning® Elplasia® plates utility for high content imaging applications, HEY-T30 and SKOV3 spheroids were exposed to varying concentrations of cisplatin for 24 hours. Spheroids were then stained with Hoechst and propidium iodide (PI) and imaged with the Thermo Scientific CellInsight™ CX7 high-content screening platform. Representative images were captured directly in Corning Elplasia plates, demonstrating ability of culture and subsequent imaging of spheroids (Figure 5). Image analysis of PI staining reveals concentration-response curves for cisplatin cytotoxicity with IC₅₀ values of 20 μM and 9.4 μM for HEY-T30 and SKOV3 spheroids respectively, which are consistent with those reported in literature (Figure 6)²⁻³.




Figure 5. High content imaging of spheroids within a Corning Elplasia plate. Representative Z-stack images of SKOV3 (top) and HEY-T30 (bottom) spheroids exposed to cisplatin (50 μ M left, 12.5 μ M middle, and 0 μ M right) for 24 hours and stained with Hoechst and propidium iodide to assess cell viability. Images taken with the Thermo Scientific CellInsight CX7 confocal imager using a 4X objective.




Figure 6. Concentration-dependent cisplatin cytotoxicity. Concentration-dependent responses of SKOV3 and HEY-T30 spheroids to 24 hours of cisplatin exposure at varying concentrations. Data shown is standard error from 3 independent studies. N=24 per dose. TC_{50} = 9.4 μ M (SKOV3) and 20 μ M (HEY-T30).

Conclusions

3D cell culture has become an important tool of the drug discovery process, which allows for research that is often more relevant than more traditional 2D cell culture. With the increase in the importance and implementation of 3D cell culture, the need for better 3D tools to generate more spheroids will also increase. Corning® Elplasia® plates provide a format to generate multiple, consistent spheroids per well for spheroid scale up or assays where more samples or larger signal is required. The black sidewalls and clear bottom design of Corning Elplasia plates make them ideally suited for imaging, luminescent, and fluorescent cell-based assays. Additionally, they are easy to use and can produce data with high signal-to-background ratios and consistent Z' values.

References

1. Imamura Yoshinori, et al. Comparison of 2D- and 3D-culture models as drug-testing platforms in breast cancer. *Oncology reports* 33(4); 1837-1843 (2015).
2. Song Jin, et al. Annexin XI is associated with cisplatin resistance and related to tumor recurrence in ovarian cancer patients. *Clinical Cancer Research* 13(22):6842-6849 (2007).
3. He Yifeng, et al. The changing 50% inhibitory concentration (IC₅₀) of cisplatin: A pilot study on the artifacts of the MTT assay and the precise measurement of density-dependent chemoresistance in ovarian cancer.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

CORNING

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

Chinese Mainland
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6572-9740
f 65 6735-2913

Taiwan

t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CSEurope@corning.com

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands

t 020 655 79 28
f 020 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400