

Corning® Matrigel® Matrix-coated Transwell® Permeable Supports for Enhancing Hepatocyte Differentiation from Human Embryonic Stem Cells

CORNING

Application Note

Miya Hsieh¹, Ronald A. Faris¹, Luke Dimasi¹, Kevin Chen²,
Lih-Tao Hsu², Hui-Ting Huang², and Wannshin Chen²

¹Corning Incorporated, Corning Life Sciences, Tewksbury, MA USA

²Biomedical Technology and Device Research Laboratories,
Industrial Technology Research Institute, Hsinchu, Taiwan

Introduction

Unforeseen adverse drug reactions in the liver is a major reason for the failure of drug candidates in clinical trials and for the after-market removal of drugs¹. These issues arise as a result of drug metabolism occurring in the body. Primary human hepatocytes are considered the gold standard model to identify potential drug-drug interactions and toxicity using *in vitro* assays; however, their limited availability and inter-individual variability constrain their routine use early in the drug development process². In recent years there has been considerable interest in generating renewable sources of hepatocytes from stem cells for drug safety studies³. Hepatic differentiation is typically accomplished by plating stem cells on tissue culture plates coated with collagen or Corning Matrigel matrix and sequentially culturing the cells with cytokines and growth factors that mimic fetal liver development⁴. While stem cell-derived hepatocytes are a promising approach to address the limitations of primary cells, existing differentiation protocols typically produce hepatic progeny whose basal and induced expression levels of major CYP drug metabolism genes (e.g., CYP1A2, CYP2B6, and CYP3A4) and key hepatic transcription factors is lower than what is observed in side-by-side comparisons with primary hepatocytes⁵. CYP3A4 is responsible for the metabolism of 40% to 50% of drugs on the market⁶. There is a clear need to identify culture formats and media that will further enhance the maturation of stem cell-derived hepatocytes.

In the liver, hepatocytes are normally arranged in one to two cell thick cords separated by sinusoids containing endothelial cells, Kupffer cells, Ito cells, and blood. Because of this unique tissue architecture, hepatocytes are polarized cells with two membrane surfaces exposed to blood⁷. We reasoned that carrying out hepatic differentiation of stem cells on permeable membrane supports, that permits exposure to nutrients, cytokines, growth factors, etc. from two directions, may enhance the maturation process. To investigate this question, hepatocyte progenitors derived from TW6 embryonic stem cells were cultured on Transwell permeable supports coated with Matrigel matrix and the basal and induced CYP expression was assayed using

Promega's P450 CellTiter-Glo® assay kits and compared to HepaRG™ cells. Our results demonstrate that the Transwell culture format enhances both the basal and pharmacologically induced expression of CYP1A2 and CYP3A4.

Material and Methods

The TW6 human pluripotent stem cells used in this study were kindly provided by Wannshin Chen (Industrial Technology Research Institute, Hsinchu, Taiwan). HepaRG and recommended medium for CYP induction assays were procured from Life Technologies. DMEM/F12 medium, RPMI medium, B-27 supplement, Insulin-selenium-transferrin, Transwell 0.4 µm PET membrane cell culture inserts (Corning Cat. No. CLS3450), Corning Hepatocyte Maintenance Medium (CHHM) (Corning Cat. No. 40-550-CV), and Corning Matrigel hESC qualified matrix (Corning Cat. No. 354277) were from Corning. Hepatocyte growth factor (HGF), fibroblast growth factor 1 (FGF1), fibroblast growth factor 4 (FGF4) were purchased from Peprotech.

Hepatic Differentiation

TW6 cells were routinely maintained on mouse embryonic fibroblast feeder layers. hESC colonies were mechanically dissected into small pieces and re-plated on 2% collagen-coated dishes. Hepatic differentiation was carried out as previously described⁸. Briefly, cells were cultured with RPMI medium containing 1X B27 and activin A (100 ng/mL) to induce endoderm differentiation followed by co-treatment with fibroblast growth factor 4 (FGF4, 10 ng/mL) and hepatocyte growth factor (HGF, 10 ng/mL) to commit endodermal precursors to the hepatocyte lineage (Figure 1). At day 8, hepatic progenitors were subcultured and seeded onto Tissue Culture (TC)-treated polystyrene plates or Transwell inserts each coated with Corning Matrigel matrix (Figure 2). The cells were further differentiated in DMEM medium containing HGF (20 ng/mL) for 5 days, followed by DMEM/F12 medium plus oncostatin M (OSM) and dexamethazone (DEX) from day 9 to day 18. At day 19, cultures were overlaid with Matrigel matrix (25 µg/well) and re-fed with CHHM containing proprietary nutrient and non-nutrient differentiation supplements for 7 days to further promote hepatocyte maturation. For CYP induction studies hepatocytes were subsequently maintained in CHHM for 3 days with or without inducers before assaying cytochrome P450 enzymatic activity.




Figure 1. Hepatic differentiation protocol

CYP Induction Studies

For CYP induction studies, TW6-derived hepatocytes were cultured in CHHM containing 100 μ M omeprazole or 10 μ M rifampicin for 72 hours to induce CYP1A2 and CYP 3A4, respectively. Control cultures were maintained in CHHM containing 0.1% DMSO. CYP1A2 and CYP3A4 enzymatic activity was measured using P450 CellTiter-Glo[®] assay kits (Promega) according to the manufacturer's instructions. The activity of TW6-derived hepatocytes was compared to HepaRG[™] cells cultured on collagen-coated TC-treated plates according to the manufacturer instructions.

Results and Discussion

The purpose of this study was to determine if differentiating hepatic progenitors on permeable supports would enhance their expression of CYP1A2 and CYP3A4. To investigate this question TW6 stem cell-derived hepatic progenitors were subcultured onto TC-treated dishes or Transwell[®] inserts coated with Corning[®] Matrigel[®] matrix on day 8. Hepatic progenitor differentiation was promoted by the sequential exposure to growth factors and cytokines which promote fetal liver development⁴. Hepatocyte maturation was further enhanced by culturing TW6 hepatocyte-like

cells in Corning hepatocyte maintenance medium supplemented with nutrients and small molecules that enhanced the expression of liver-related transcription factors and decreased the expression of alpha-fetoprotein in 2D culture (data not shown).

In this study we observed that differentiating hepatic progenitors on Transwell supports in CHHM plus supplements further enhanced their maturation as measured by the enhanced expression of two key drug metabolism genes, CYP1A2 and CYP3A4 (Figures 3 and 4). Importantly, both the basal and induced expression of both CYP enzymes were elevated 2- to 3-fold in comparison to cells differentiated on Matrigel matrix-coated TC-treated plates. Moreover, CYP function of hepatocytes generated on Transwell inserts was similar to the activity of HepaRG cells maintained according to supplier's recommendations.

The results from this study clearly indicate that differentiating hepatic progenitors on permeable membrane supports enhances their subsequent maturation. How the Transwell culture format promotes hepatocyte differentiation is presently unclear. Other investigators have shown that co-culturing hepatocytes and endothelial cells on opposite sides of the Transwell membrane prolongs cell viability and hepatic function. The authors propose that their co-culture model mimics the *in vivo* liver sinusoid




Figure 2. Differentiation of hepatic progenitors on Matrigel-coated Transwell inserts




Figure 3. Differentiation of hepatic progenitors on Corning Matrigel matrix-coated Transwell inserts enhanced both basal and rifampin-induced CYP3A4 expression.





Figure 4. Differentiation of hepatic progenitors on Corning Matrigel matrix-coated Transwell inserts enhanced both basal and omeprazole-induced CYP1A2 expression.

architecture⁹. More recently, Dash, et al. reported that Transwell culture in combination with hemodynamic flow minimizes the rapid dedifferentiation of primary cultured rat hepatocytes and improves the expression of Phase I and Phase II drug metabolism genes¹⁰. In our study differentiated hepatocytes form a confluent monolayer on the Corning® Matrigel® matrix-coated Transwell® inserts, thus enabling exposure to nutrients and oxygen from two directions (Figure 2). This culture format is readily amenable for drug screening assays.

References

- Shuster D, Laggner C, and Langer T. Why drugs fail—A study on side effects in new chemical entities. *Curr. Pharm. Des.* 11:3545-3559 (2005).
- eCluysse EL. Human hepatocyte culture systems for in vitro evaluation of cytochrome P450 expression and regulation. *Eur. J., Pharm. Sci.*, 13:343-368 (2001).
- Szkolnicka D, Zhou W, Lucendo-Villarin B, and Hay DC. Pluripotent Stem Cell-Derived Hepatocytes: Potential and Challenges in Pharmacology. *Annu. Rev. Pharmacol. and Toxicol.* 53:147-159 (2013).
- Schwartz RE, Fleming HE, Khetani SR, and Bhatia SN. Pluripotent stem cell-derived hepatocyte-like cells. *Biotechnol. Adv.* 32:504-513 (2014).
- Godoy P, Schmidt-Heck W, Natarajan K, Lucendo-Villarin B, Szkolnicka D, Asplund A, Björquist P, Widera A, Stöber R, Campos G, Hammad S, Agapios Sachinidis, Chaudhari U, Damm G, Weiss TS, Nüssler A, Synnergren A, Edlund K, Küppers-Munther B, Hay DC, and Hengstler JG. Gene networks and transcription factor motifs defining the differentiation of stem cells into hepatocyte-like cells. *J Hepatol.* 63:934-942 (2015).
- Zhou S, Yung Chan S, Cher Goh B, Duan W, Huang M, McLeod HL. Mechanism-based inhibition of cytochrome P450 3A4 by therapeutic drugs. *Clin. Pharmacokinet.* 44:279-304 (2005).
- Hardman RC, Volz DC, Kullman SW, and Hinton DE. An in vivo look at vertebrate liver architecture: Three-dimensional reconstructions from Medaka (*Oryzias latipes*). *The Anatomical Record* 290:770-782 (2007).
- Chen-Ming Chen, Wannhsin Chen, and Lih-Tao Hsu. Cultured hepatocyte and method for preparing the same. US Patent Application 2015/0166952 (June 18, 2015).
- Kang YB (Abraham), Rawat S, Cirillo J, Bouchard M, Noh H (Moses). Layered long-term co-culture of hepatocytes and endothelial cells on a Transwell membrane: Toward engineering the liver sinusoid. *Biofabrication.* 2013;5(4):045008. doi:10.1088/1758-5082/5/4/045008.
- Dash A, Simmers MB, Deering TG, Berry DJ, Feaver RE, Hastings NE, Pruett TL, LeCluysse EL, Blackman BR, and Wamhoff BR. Hemodynamic flow improves rat hepatocyte morphology, function and metabolic activity in vitro. *Am. J. Cell Physiol.* 304:C1053-C1063 (2013).



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Customers outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476
www.corning.com/lifesciences

Worldwide Support Offices

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6572-9740
f 65 6861-2913

Taiwan

t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands

t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

All Other European Countries

t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA grupoLA@corning.com

Brasil
t (55-11) 3089-7400
Mexico
t (52-81) 8158-8400

CORNING

FALCON

AXYGEN

GOSSSELIN

PYREX

For a listing of trademarks, visit www.corning.com/clstrademarks.
All other trademarks are the property of their respective owners.