

Drop Test Performance of Corning® PET and Competitor PETG Bottles

Application Note

CORNING

*Ludivine Rineau
Corning Incorporated, Life Sciences
CETC Fontainebleau, France
METROPACK, Reims, France*

Introduction

Polyethylene terephthalate (PET) and Polyethylene Terephthalate Glycol (PETG) bottles are both recommended for storage of media and sera for cell culture applications. This study sought to determine if there is any difference in the mechanical performance of PET and PETG bottles when used to store frozen solutions down to -70°C.

Materials and Methods

Bottle Samples

- ▶ Corning PET octagonal bottle, 1L (Corning Cat. No. 432334)
- ▶ Competitor PETG square bottle, 1L

Sample Preparation

Thirty samples of each bottle type were filled with 1L of solution comprised of 1 g/L of glucose and 9 g/L of NaCl. This solution was used to mimic typical cell culture preparations. The bottles were capped according to the manufacturer's torqueing specifications using a Steinfurth torque measuring instrument. The bottles were placed upright in a -70°C freezer for 48 hours and then transferred to a -40°C freezer for 48 hours.

Drop Testing

The bottles were removed from the -40°C freezer and immediately dropped one time from 910 mm (35.8 inches) onto a steel surface attached to a concrete block in accordance with ISO 2248 and ASTM D 5276 standards. Bottles were dropped in a controlled orientation, with a large side face making initial impact, using a Lansmont drop tester. After the drop test, a visual inspection of the bottle was performed to assess bottle integrity while still frozen and after thawing.

Results

No failure or leakage was observed after drop testing frozen Corning PET octagonal bottles or the competitor PETG square bottles.

Summary and Conclusions

The mechanical performance of Corning PET octagonal bottles is comparable to that of the competitor PETG square bottles for storing frozen liquid. Corning PET octagonal bottles are suitable for storage of frozen media, sera, and other aqueous solutions down to -70°C.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea

t 82 2-796-9500
f 82 2-796-9300

Singapore

t 65 6572-9740
f 65 6861-2913

Taiwan

t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CSEurope@corning.com

France

t 0800 916 882
f 0800 918 636

Germany

t 0800 101 1153
f 0800 101 2427

The Netherlands

t 020 655 79 28
f 020 659 76 73

United Kingdom

t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA

grupoLA@corning.com

Brasil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400