

Counting Murine CAR T Cells Using the Corning® Cell Counter

Application Note

CORNING

Jeff Partridge
Corning Incorporated, Life Sciences
Bedford, MA USA

Introduction

About CAR T cells

Chimeric Antigen Receptor, or CAR T cells, are genetically modified autologous T cells that are currently used for immunotherapy to treat certain cancers¹⁻³. These T cells are modified to express CARs, which are proteins that allow T cells to recognize an antigen on targeted tumor cells⁴.

Early outcomes from CAR T-cell trials have generated impressive results in patients with refractory or relapsed leukemias and lymphomas¹⁻², and research is underway on using CAR T cells in solid tumors^{2,4}.

Generating CAR T cells

Whole blood is collected from a patient (in this application note, mouse) and PBMC (peripheral blood mononuclear cells) are separated and saved while the remaining blood cells and components are returned to the body (leukapheresis). Next, T cells are enriched from the PBMC population using various methods (e.g., density gradient, size-based fractionation, antibody bead conjugates). The T cells are genetically modified with genes of interest, such that they now express CARs, allowing them to target and destroy tumor cells. CAR T cells are activated and expanded (*ex vivo*), and in production scale there are additional QA and transportation/delivery steps. Finally, CAR T cells are infused back to the patient³. This may be indirectly to the peripheral circulation or directly to the affected site depending on the tumor type⁵.

Counting of CAR T cells

CAR T cells are often counted using time consuming manual methods such as a hemocytometer. While the hemocytometer has been considered the gold standard for cell counting since the 18th century⁶, results are frequently subjective, variable, and small differences in cell diameters cannot easily or reproducibly be discerned.

The Corning Cell Counter offers an alternative method to quickly and accurately count these cells. The counting algorithm was developed to identify and remove debris and account for dead cells, while the functionality of the application allows for the population to be gated to remove unmodified T cells from the final count.

Materials and Methods

Murine CAR T cells obtained just prior to infusion were counted on the Corning Cell Counter, then immediately manually counted and viability determined using a hemocytometer and the trypan blue exclusion method. Representative images of the area counted (3 mm²) are shown in Figure 1.

The experimental cell suspensions counted were typically heterogeneous samples, each having many dead cells, debris, and non-modified T cells, which are smaller than CAR T cells⁷.




Figure 1. Representative the Corning Cell Counter images after focusing and after counting. (A) Properly focused image before count. (B) Results of count. Live cells are circled green; dead cells are circled red. Small cells and debris may be excluded from counting by the algorithm and by user specified gating (example in Figure 2). Note that a counting chamber having counting grids may be used on the Corning Cell Counter, but grids are not necessary. (C) Magnified view of B panel. Not all cells on screen were counted. **This is normal for this application.** Many of these cells were not modified (i.e., are not CAR-T cells) and should be excluded by size and/or cell morphology⁷.

Results and Discussion

Using the *post hoc* histogram capability in the Axion app or Axion cloud, the bottom gate was set at 8 μm (Figure 2) due to the size difference between unmodified T cells and CAR T cells. For the experimental data generated during this evaluation, the results over nine independent data sets are shown in Figures 3 and 4.

When comparing the two counting methods using a Paired t-test, p values for both cell densities and percent viability were greater than 0.05, indicating there were no statistical differences between each data set. These results indicate the suitability of the Corning® Cell Counter for accurate counting of CAR T cells.

Although the human eye is quite sensitive, it may not be able to tell the difference between cells of five to eight microns, while the Corning Cell Counter can. Thus, gating out based on known cell size is a simple process with the Corning Cell Counter but not as easy when manually counting. This functionality allowed for accurate counts of murine CAR T cells by a trained operator.

Conclusions

- ▶ The Corning Cell Counter may be used in murine CAR T production where accurate cell counting is required. The Corning Cell Counter is intended as a replacement method for manual counting with a hemocytometer.
- ▶ While the Corning Cell Counter is very easy to use, proper focusing techniques need to be observed to avoid subjectivity in counts from different users. Increased familiarity with the Corning Cell Counter improves this accuracy. With every count, users will gain experience and proficiency.
- ▶ It is possible to count cell suspensions of a higher density on the Corning Cell Counter than by manual means. Thus the overall cell population will more accurately be represented, and the associated problem of “clicker fatigue” is avoided altogether.

References

1. Wang X, et al. Clinical Manufacturing of CAR T Cells: Foundation of a Promising Therapy. *Molecular Therapy — Oncolytics* (2016) 3, 16015.
2. Levine BL, et al. Global Manufacturing of CAR T Cell Therapy. *Molecular Therapy: Methods & Clinical Development* Vol. 4, March 2017:92-101.
3. Fesnak AD, et al. Production of Chimeric Antigen Receptor T cells. *Nature Protocols*, Document # 27041, Version 1.0.0 (2017). Accessed online from www.nature.com/nprot/posters/cart/index.html
4. <https://www.lls.org/treatment/types-of-treatment/immunotherapy/chimeric-antigen-receptor-car-t-cell-therapy>
5. <https://www.mskcc.org/announcements/combination-immunotherapy-boosts-car-cell-approach-solid-tumors>
6. Lin-Gibson S, et al. Summary of the National Institute of Standards and Technology and US Food And Drug Administration cell counting workshop: Sharing practices in cell counting measurements *Cytotherapy*, 2018; 20:785-795.
7. Personal communication.




Figure 2. Histogram showing bottom gate set at 8 μm cell diameter. Cells below this threshold were excluded from the count. Dead cells (red bars) and live cells (green bars) can be toggled to independently examine each subpopulation.




Figure 3. Cell densities by count method. Paired t-test $p = 0.135$.




Figure 4. Cell viabilities by count method. Paired t-test $p = 0.274$.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA

t 800.492.1110
t 978.442.2200

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

LATIN AMERICA

grupoLA@corning.com

Brazil

t 55 (11) 3089-7400

Mexico

t (52-81) 8158-8400