

Corning® 100 mm Aseptic Transfer Cap for 5L Erlenmeyer Flasks Provides an Aseptic Solution for Liquid Transfer and Can Be Utilized Throughout Cell Culture

CORNING

Application Note

Hilary A. Sherman, Lori Romeo, and Hannah Gitschier
Corning Incorporated, Life Sciences
Kennebunk, ME USA

Introduction

Closed cell culture systems play a significant role in reducing the risk of contamination during liquid handling. The Corning 100 mm aseptic transfer cap (ATC) offers a closed system solution for transferring liquid aseptically into or out of a Corning 5L Erlenmeyer flask. The ATC also contains a hydrophobic vent with 0.2 µm pore size for gas exchange, making the cap ideal for use during cell culture. The flexible tubing has a male medical plastic coupler (MPC) and is also heat sealable, allowing for a multitude of connectivity options. The vessel can be filled via gravity flow, peristaltic pump, or with negative pressure applied to the 0.2 µm vent filter of the ATC. To empty the vessel, a peristaltic pump can be used, or positive pressure can be applied to the transfer cap vent filter. Here we demonstrate the use of the ATC to seed and expand two suspension cell lines, 5/9 m alpha3-18 and Sf9 cells in the 5L Erlenmeyer flasks.

Materials and Methods

5/9 m alpha3-18 cells (ATCC® Cat. No. CRL-10154) were cultured in CD OptiCHO™ (Thermo Fisher Cat. No. 12681-011) supplemented with 8 mM L-Glutamine (Corning Cat. No. 25-005-CI) and 1X Hypoxanthine and Thymidine (Corning Cat. No. 25-047-CI). Cells were incubated in an Infors HT Multitron Pro orbital shaker at 37°C and 5% CO₂. 5/9 m alpha3-18 cells were seeded into Corning 5L Erlenmeyer flasks (Corning Cat. No. 431685) and 5L Erlenmeyer flasks with 100 mm aseptic transfer caps (Corning Cat. No. 11500) at a density of 3 x 10⁵ cells/mL and a volume of 2.5 L per vessel. Cells and medium were poured directly into the 5L Erlenmeyer flasks without ATCs, and transferred into the 5L Erlenmeyer flasks with ATCs utilizing the ATC by way of a peristaltic pump at a flow rate of approximately 1L/min. Once seeded, cells were rotated at 90 rpm and counted daily, for 4 days, utilizing the Beckman Coulter Vi-CELL™ Cell Viability Analyzer. The study was repeated one additional time for a total of 5 vessels with ATCs and 6 vessels with standard caps.

Sf9 cells (Thermo Fisher Cat. No. 12659-017) were cultured in Sf-900™ II SFM (Thermo Fisher Cat. No. 10902104) at a density between 1 x 10⁶ to 8 x 10⁶ cells/mL in a Gallenkamp floor model orbital shaker set to 28°C. Sf9 cells were seeded into Corning 5L Erlenmeyer flasks and 5L Erlenmeyer flasks with 100 mm aseptic transfer caps at a density of 1 x 10⁶ cells/mL in 3.5 L per vessel. Cells and medium were poured directly into the 5L Erlenmeyer flasks without ATCs or transferred utilizing the ATC by applying a vacuum of 2.5 inches of mercury to the ATC vent filter. Cells were rotated at 90 rpm and counted daily on the Nova BioProfile® FLEX analyzer for 4 days. The study was repeated one additional time for a total of 6 vessels per condition.

Results and Discussion

5L Erlenmeyer flasks affixed with 100 mm aseptic transfer caps allow for easy fluid transfer into and out of the vessels while also reducing the risks for contamination typically acquired by pouring. 5/9 m alpha3-18 cells seeded via an ATC at a 1L/min. flow rate from a peristaltic pump, and cultured with the ATC in place, showed equivalent high viability and cell growth throughout the 4-day culture period when compared to standard cap vessels (Figure 1). Sf9 cells seeded via an ATC with 2.5 inches of mercury negative pressure applied to the ATC vent filter, and cultured with the ATC in place, also displayed equivalent high viability and cell growth compared to standard cap vessels (Figure 2). Our data demonstrate that fluid can be transferred via peristaltic pump or by application of negative pressure without a detrimental effect to the culture, and that the cells thrive with high viability with the ATC left in place throughout the duration of culture.



Figure 1. 5/9 m alpha3-18 cell viability was maintained above 96% regardless of whether a Corning aseptic transfer cap (ATC) or standard cap was used (right axis). Throughout the four-day culture there was no statistical difference in cell densities achieved between vessels with or without an ATC (left axis). Two-way ANOVA $p > 0.05$, $n = 5$ for ATC vessels and 6 for standard cap vessels.




Figure 2. Sf9 cell viability was maintained above 97% regardless of whether a Corning aseptic transfer cap (ATC) or standard cap was used (right axis). Throughout the four-day culture there was no statistical difference in cell densities achieved between vessels with or without an ATC (left axis). Two-way ANOVA $p > 0.05$, $n = 6$ for all vessels.

Conclusions

- ▶ Corning 100 mm aseptic transfer caps can be utilized to aseptically transfer liquid and cell cultures into Corning 5L Erlenmeyer flasks by way of peristaltic pump or vacuum pump.
- ▶ Utilization of the ATC for filling 5L Erlenmeyer flasks has no detrimental effect on cell viability or cell growth when compared to standard cap vessels that were seeded by pouring.
- ▶ There was no statistical difference in cell densities or viability from cells cultured in 5L Erlenmeyer flasks with ATCs compared to standard 100 mm cap vessels.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6735-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE
CSEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t 55 (11) 3089-7400

Mexico
t (52-81) 8158-8400