

# Corning® Lambda™ EliteMax Benchtop Pipettor for ELISA Workflows

CORNING

## Application Note

Hilary Sherman and Ann Ferrie  
Corning Incorporated, Life Sciences  
Kennebunk, ME USA

### Introduction

Enzyme-linked immunosorbent assays or ELISAs are an essential analytical tool for detecting and quantifying molecules in solution<sup>1</sup> which is why they are frequently used in drug discovery, diagnostic tests, and quality control applications<sup>2</sup>. Running an ELISA requires multiple and tedious liquid handling steps. From adding proteins or antibodies to multiple washes to remove unbound components to the addition of the detection agent requires a tremendous amount of liquid manipulations. Here we demonstrate use of the Corning Lambda EliteMax benchtop automated pipettor to easily perform the numerous liquid handling steps needed for an indirect ELISA.


### Materials and Methods

All liquid addition and removal steps were performed using plate-to-plate transfer protocols at medium speed and Axygen® 96-well tips (Corning FX-250-L-R). 150 µL of serially diluted recombinant human insulin (MP Biomedicals 193900), with a starting concentration of 0.4 µg/mL, from an 8-channel trough (Corning RES-MW8-HP-SI) were added across each well of 3 Corning 96-well clear flat bottom polystyrene high bind microplates (Corning 9018). Microplates were then sealed using Corning 96-well microplate aluminum sealing tape (Corning 6570) and stored at 2-8°C overnight prior to use. When ready, plates were removed from 2-8°C and allowed to warm to room temperature. Meanwhile, a 1X PBS wash buffer, containing 0.1% TWEEN® 20 (Sigma P2287), was prepared from a 10X phosphate buffered saline (PBS; Corning 46-013-CM) stock solution diluted with Type II water. The serially diluted insulin was removed and replaced with 150 µL of blocking buffer consisting of 3% bovine serum albumin (Sigma A9576-50ML) in 1X PBS wash buffer and incubated for 1 hour at room temperature. After the blocking buffer was removed, plates were washed 3 times with 200 µL of 1X PBS wash buffer by linking a plate-to-plate transfer protocol to remove liquid and a plate-to-plate transfer protocol to add 1X wash buffer. For ease of washing, Axygen single well reagent reservoirs with 96-bottom troughs (Corning RES-SW96-LP) were used for the 1X PBS wash buffer as well as for waste. Next, 100 µL of 0.76 µg/mL Mouse IgG1 insulin monoclonal antibody (ThermoFisher MA1-83252), diluted in 1X PBS wash buffer, was added to each well and incubated at room temperature for 1 hour. Plates were then washed three additional times as previously described. After washing the wells of unbound primary antibody, 150 µL of 0.5 µg/mL goat anti-mouse IgG (H+L) poly-HRP secondary antibody, HRP (ThermoFisher 32230) diluted in 1X PBS wash buffer was added to each well and incubated at room temperature for 1 hour. Microplates were washed 3 times after

secondary antibody incubation following the previously described method. Finally, 100 µL of 1:1 mixed TMB 2-Component Microwell Peroxidase Substrate kit (SeraCare 5120-0047) was added to each well and incubated for approximately 5 minutes at room temperature. This color changing reaction was stopped with the addition of 100 µL per well of Stop solution (TMB 5150-0021) as per manufacturer's protocol. The change in color in the wells were then read on the Envision® multimode plate reader (PerkinElmer 2105-0010) at 450 nm.

### Results and Discussion

Use of automation for repetitive processes can speed up assay workflows and reduce ergonomic risk but assay robustness and repeatability must not be compromised during implementation. Three independent ELISA assays were performed (n = 12) fully utilizing the Lambda EliteMax for all liquid handling tasks resulting in low intra-plate and inter-plate variability (Figure 1).



**Figure 1. Dose dependent insulin binding achieved from ELISA automated with Corning Lambda EliteMax benchtop pipettor.** Dose-dependent insulin binding curves generated from 3 microplates shown with mean and standard error. N = 12 per concentration. Average EC<sub>50</sub> from 3 studies is 0.025.

### Conclusions

ELISAs can be very labor-intensive and time-consuming due to numerous liquid addition and removal steps. Automating this process can save time, reduce intra-plate variability, and improve end user fatigue by reducing manual pipetting and improving ergonomics. Here, we demonstrated use of the Corning Lambda EliteMax benchtop pipettor to fully automate liquid handling steps for an indirect ELISA.

## References

1. Aydin S. A short history, principles, and types of ELISA, and our laboratory experience with peptide/protein analyses using ELISA. *Peptides* 72 (2015): 4-15.
2. Hosseini S. et al. General overviews on applications of ELISA. *Enzyme-linked Immunosorbent Assay (ELISA)*. Springer, Singapore, 2018. 19-29.

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use or general laboratory use only.\* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. \*For a listing of US medical devices, regulatory classifications or specific information on claims, visit [www.corning.com/resources](http://www.corning.com/resources).

*Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.*

# CORNING

Corning Incorporated  
Life Sciences

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

## NORTH AMERICA

t 800.492.1110  
t 978.442.2200

## ASIA/PACIFIC

### Australia/New Zealand

t 61 427286832

### Chinese Mainland

t 86 21 3338 4338

## India

t 91 124 4604000

## Japan

t 81 3-3586 1996

## Korea

t 82 2-796-9500

## Singapore

t 65 6572-9740

## Taiwan

t 886 2-2716-0338

## EUROPE

CSEurope@corning.com

## France

t 0800 916 882

## Germany

t 0800 101 1153

## The Netherlands

t 020 655 79 28

## United Kingdom

t 0800 376 8660

## All Other European Countries

t +31 (0) 206 59 60 51

## LATIN AMERICA

grupoLA@corning.com

## Brazil

t 55 (11) 3089-7400

## Mexico

t (52-81) 8158-8400