Stem Cell Therapy Production

Seeding, Expanding, and Harvesting Stem Cells

Key considerations for working with MSCs, iPSCs, and NSCs through different stages of the workflow

Mesenchymal Stem Cells (MSCs)

Growth as Individual Cells Vessel selection

- Traditional: Dishes, plates, T-flasks
- Stacked vessels: Corning CellSTACK, and HYPERStack vessels
- Bioreactor: Corning Ascent Fixed Bed Reactor (FBR), microcarriers in bioreactor

Growth surface selection

- Surface treatment
 - TC-treated
 - Corning CellBIND
- Surface coating (pre-coated or self-coated)
 - Collagen
 - Fibronectin
 - Corning Synthemax II

Seeding Density 200-12,000 cells/cm² (most commonly 1,000-6,000 cells/cm²)

Lower seeding density

- Increased proliferation potential/fold-expansion
- Fewer passages to reach target yield

Higher seeding density

- Reduced time to reach target cell density
- Economical for low output
- Increased stress to cells due to paracrine signaling, leading to stress fibers

Doubling Time 24-40 hours

Passaging Time 3-7 days

Target Confluency 75-80%

Media Change Every 2-3 days





Culture options and neural differentiation paths

Induced Pluripotent Stem Cells (iPSCs)

Growth as Clusters

Substrate on dishes, plates, T-flasks, or CellSTACK vessels

Mouse embryonic fibroblasts (MEFs)

- Irradiated animal cells - Safety concerns

Corning Matrigel matrix

- From mouse sarcoma cells
- Not fully defined

Corning Synthemax II vitronectin substrate

- Synthetic
- Xeno-free

Seeding Density 10,000-20,000 cells/cm²

Doubling Time 16-20 hours

Passaging Time 4-5 days

Passaging Criteria


60-75% confluency and/or medium-to-large colony size and/or signs of spontaneous differentiation

Microscopic observation daily for

- iPSC-like morphology
- Differentiated cells







2D Coculture with stromal cells

Corning rLaminin-521 (human)

- Recombinant human Laminin protein

Growth as Individual Cells

Recommend with Ascent FBR or microcarriers and bioreactor



- Confluency

Accutase cell detachment solution


- Gentle, enzyme-free dissociation preserves genomic stability

Corning CellStripper solution

- Non-enzymatic cell dissociation solution formulated with a proprietary mixture of chelators

Manual passaging with pipet tip or cell scraper

Media change Daily



Neural Stem Cells (NSCs)

Growth at High Density

Growth surface to promote attachment on dishes, plates, T-flasks, or CellSTACK vessels

- TC-treated
- CellBIND


Add positive charge

- Poly-L-Ornithine
- Poly-L-Lysine
- Poly-D-Lysine
- Corning PureCoat Amine
- Add Extracellular Matrix (ECM)
- Laminin

Seeding Density 20,000-100,000 cells/cm²

Doubling Time 20-48 hours (very limited before differentiation)

Confluency at Passaging 95-100%



Neurospheres

3D

3D Brain organoids

Learn more about the most important considerations for working with different stem cell types — MSCs, iPSCs and NSCs — through different stages of the workflow.

www.corning.com/celltherapy

CORNING

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning, Accutase, Ascent, CellBIND, CellStack, CellStripper, HYPERStack, Matrigel, PureCoat, and Synthemax are trademarks of Corning Incorporated. For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners. ©2023 Corning Incorporated. All rights reserved. 9/23 CLS-AN-763