Assessing the particle concentration of mRNA-LNP using Corning[®] Videodrop

A new metric for Lipid Nanoparticles (LNP) manufacturing, optimization and standardization

Application Note

Ivan Ciganek¹, Christophe Delehedde¹, Claire Counil¹, Thomas Ador¹, Anthony Delalande¹, Chantal Pichon¹, and Marie Berger² ¹CBM and University of Orléans and CNRS; ²Myriade

Introduction

Corning Videodrop is an innovative nanoscale imaging technology. Based on the principle of interferometry, Corning Videodrop makes it possible to measure the size and concentration of nanoparticles.

- In real-time (40 s)
- In a single drop (7 μL)
- Between 80 to 500 nm
- ▶ Within a concentration range of 10⁸ to 10¹⁰ part/mL
- Without labeling

Here, we propose to use Corning Videodrop for the analysis of mRNA-LNP and highlight the interest of measuring the particle concentration of mRNA-LNP (Figure 1).

The current analytical strategies employed for mRNA-LNP characterization lie primarily on size, polydispersity, and zetapotential measurements. However, particle concentration, evaluated through a single particle measurement technique, is an important parameter to monitor for several reasons (Figure 2).

- Quality control tests
- Stability assessment
- Standardization for comparative studies

When developing new formulations, manufacturing processes, or different storage conditions, quantifying the number of particles enables efficient comparison studies. By comparing particle numbers under different conditions, researchers can assess the impact of the process or composition on transfection activity.



Figure 1. Importance of LNP particle concentration.



Figure 2. Corning Videodrop applications.



1. Solution Preparation RNA/Lipid – Lipid Composition Optimization

Two LNPs with different lipid compositions were produced using the same process. Particle concentration was not affected by the lipid composition (Figure 3); however, Lipid Composition 1 exhibits higher activity compared to Lipid Composition 2.

Despite having approximately the same number of particles, Lipid Composition 1 was more efficient compared to Lipid Composition 2.

2. LNP Formulation – Mixing Speed Optimization

We compared LNPs having the same composition but prepared with different mixing speeds (flow rates) on the microfluidic device.

The mixing speed step had a significant influence on the concentration of the generated particles (Figure 4).

The activity was directly related to the number of particles.

3. Buffer Exchange and Concentration – Buffer Exchange Optimization

Buffer exchange and concentration steps are crucial and require optimization to improve the yield.

This experiment examined the influence of incorporating a buffer exchange step.

The buffer exchange step resulted in a 75% reduction in the number of particles (Figure 5).

4. Fill and Finish – Storage Optimization

To ensure product quality, it is essential to optimize storage conditions.

We tested the same product stored at +4°C or -20°C.

The results indicated that the storage at -20°C led to a decrease in the number of particles as well as a loss in activity (Figure 6).



Figure 3. Comparision of the two LNPs with different lipid compositions.

Composition 1 Composition 2



Figure 4. Influence of the mixing speed on the concentration of the generated particles.



Figure 5. Influence of incorporating a buffer exchange step.



Figure 6. Influence of the storage temperature on the number of particles.

Number of particles*

4

0

Activity

*Number of particles is calculated by multiplying the measured concentration by the volume used for activity assays.

1200

0

Conclusions

- Corning[®] Videodrop is a ready-to-use tool for measuring the particle concentration of LNPs in a single drop.
- Corning Videodrop enables better control of the production process by quickly and easily measuring the concentration.
- Corning Videodrop is the ideal tool to complete the analytical strategy to assess LNP quality.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated Life Sciences

www.corning.com/lifesciences

NORTH AMERICA t 800.492.1110 t 978.442.2200

ASIA/PACIFIC Australia/New Zealand t 61 427286832 Chinese Mainland t 86 21 3338 4338 India t 91 124 4604000 Japan t 81 3-3586 1996 Korea t 82 2-796-9500 Singapore t 65 6572-9740 Taiwan t 886 2-2716-0338 EUROPE CSEurope@corning.com France t 0800 916 882 Germany t 0800 101 1153 The Netherlands t 020 655 79 28 United Kingdom t 0800 376 8660 **All Other European Countries** t +31 (0) 206 59 60 51

LATIN AMERICA grupoLA@corning.com Brazil t 55 (11) 3089-7400 Mexico t (52-81) 8158-8400