

Corning® MSCulture Max™ Media



Instructions for Use

Media Basics

Corning MSCulture Max media requires a basal medium and supplement that have been optimized for expansion of mesenchymal stem/stromal cells (MSCs). Both components must be combined to make MSCulture Max media complete as outlined in the Preparation Method section.

Ordering Information

Corning MSCulture Max Media

NOTE: To achieve optimal MSC expansion performance, you will need to order both the basal medium and the corresponding supplement individually.*

Cat. No.	Description	Shelf Life	Shelf Life of Mixed Basal Medium and Corresponding Supplement*
Corning MSCulture Max Media, Standard Research-grade			
42-010-CV	MSCulture Max basal medium	18 months	2 weeks
42-000-CR	MSCulture Max supplement	2 years	
Corning MSCulture Max-XF Media, Xeno-free Research-grade			
42-010-CV	MSCulture Max basal medium	18 months	2 weeks
42-100-CR	MSCulture Max-XF supplement	2 years	
Corning MSCulture Max-GMP Media, Xeno-free GMP-grade for Bioproduction			
42-010-CV	MSCulture Max basal medium	18 months	1 week
42-200-CR	MSCulture Max-GMP supplement	3 years	

*As stated above, both components must be combined to make MSCulture Max media complete per the Preparation Method section below to achieve optimal MSC expansion performance.

Storage and Shipping Information

Basal medium should be stored at 2°C to 8°C (protected from light), while the supplement should be stored at -35°C to -15°C. Components are shipped separately based on temperature requirements and may arrive in two different shipments on different dates.

Preparation Method

1. Bring basal medium and corresponding supplement to room temperature (protected from light).
2. Prepare complete media by aseptically combining basal medium with corresponding supplement.
3. Mix well, and use immediately or store at 2°C to 8°C per shelf life listed above.

How to Purchase: For specific availability in your region and purchasing options, terms and conditions of sale, customer/product support, and certificates, visit www.corning.com/how-to-buy.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand

t 61 427286832

India

t 91 124 4604000

Japan

t 81 3-3586 1996

Korea

t 82 2-796-9500

Singapore

t 65 6572-9740

Taiwan

t 886 2-2716-0338

EUROPE

CEurope@corning.com

France

t 0800 916 882

Germany

t 0800 101 1153

The Netherlands

t 020 655 79 28

United Kingdom

t 0800 376 8660

All Other European Countries

t +31 (0) 206 59 60 51

The information contained within is accurate as of the date of publication and subject to change without notice.

For a listing of trademarks, visit www.corning.com/trademarks. All other trademarks are the property of their respective owners.

© 2026 Corning Incorporated. All rights reserved. 2/26 CLS-AN-899DOC REV1 ROW