

Corning® High Content Screening Microplates

Frequently Asked Questions

CORNING

1. I find some areas within wells where there are no cells attached. How can I avoid this?

Use of a manual pipettor can sometimes result in the formation of air bubbles which will prevent cells from reaching the microplate surface. Use an automated pipettor or liquid handler to reduce the chances of pipetting difficulties. You can also pre-wet the wells with small volumes of medium (25 µL per well for 96-well half area (A/2) microplates or 10 µL per well for 384-well microplates) and quickly centrifuge to ensure media reaches the bottom of all the wells. Cells can be manually seeded more easily immediately following this.

2. I noticed an uneven distribution of cells in the outer wells of my microplate. How do I reduce this edge effect?

After seeding the microplates, let them sit for 30 to 45 minutes in order for cells to settle prior to moving to the incubator. This will negate any effect from temperature differences from the edge of the microplate compared to the center of the microplate that occurs upon initial incubation.

3. Why are the wells round?

The round well geometry prevents wicking of liquid within the well and prevents cells from aggregating into square well corners.

4. What are the main dimensions of the microplate that I need in order to use Corning high content screening microplates with my in-house imager?

Microplate Dimensions

Representative Microplate	Cat. No.	Well Depth (mm)	Well Diameter (Top/Bottom) (mm)	Plate L x W x H (mm)	A1 Row Offset (mm)	A1 Column Offset (mm)	Well Center to Well Center Spacing (mm)	Flange or Skirt Height (mm)	Stack Height (mm)	Well Bottom Elevation (mm)	Well Bottom Thickness (mm)	Well Bottom Area (cm ²)	Distance to Bottom of Microplate (mm)
96-well A/2 film bottom	4680	8.94	5.00/4.49	127.8 x 85.5 x 14.2	14.4	11.24	9	6.096	12.49	2.4	0.127	0.159	2.349
96-well A/2 glass bottom	4580	8.94	5.00/4.49	127.8 x 85.5 x 14.2	14.4	11.24	9	6.096	12.49	2.4	0.2	0.159	2.349
384-well film bottom	4681	10.92	3.4/2.8	127.8 x 85.5 x 14.2	8.99	12.13	4.5	6.096	13.15	3.3	0.127	0.0615	1.5
384-well glass bottom	4581	10.92	3.4/2.8	127.8 x 85.5 x 14.2	8.99	12.13	4.5	6.096	13.15	3.3	0.2	0.0615	1.5

5. Why do the Corning high content screening microplates have smaller wells than other manufactures?

Corning 96-well half area (A/2) and 384-well microplates have smaller wells to allow for less use of reagents and cells, which may result in a cost savings. The 96- and 384-well bottom surface area is sufficient to capture quality data for high content imaging. Using a Thermo Scientific CellInsight™ imager, we determined the number of fields of view that could be obtained without capturing the edge of the well. For most users, these fields of view are acceptable to capture more than 1,000 nuclei per well.

Microplate Format	Fields of View	
	10X Objective	20X Objective
96-well half area (A/2)	15	121
384-well	9	36

6. I seem to have greater evaporation along the edges of my microplate, and I worry that it could affect my assay. How can I reduce this type of edge effect?

Make sure there is a sufficient media volume per well. 40 µL per 384-well microplate and 100 µL per 96-well half area (A/2) microplate should have limited evaporation for a 24- to 48-hour culture period. For longer culture periods, breathable tape (Corning Cat. No. 3345) can significantly reduce evaporation without inhibiting gas exchange.



www.corning.com/lifesciences/solutions

At Corning, we continuously strive towards improving efficiencies and developing new products and technologies for life science researchers. From assay preparation to storage, our technical experts understand your challenges and your increased need for high-quality products.

It is this expertise, plus a 160-year legacy of Corning innovation and manufacturing excellence, that puts us in a unique position to be able to offer a beginning-to-end portfolio of high-quality, reliable life sciences consumables.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, email ScientificSupport@corning.com, visit www.corning.com/lifesciences, or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

**Corning Incorporated
Life Sciences**
836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476
www.corning.com/lifesciences

**Worldwide
Support Offices**

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

**All Other European
Countries**
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA
grupoLA@corning.com
Brasil
t (55-11) 3089-7400
Mexico
t (52-81) 8158-8400

CORNING | **FALCON** | **AXYGEN** | **GOSELIN** | **PYREX**

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.