

# High Throughput Genomic and Proteomic Sample Preparation Using Corning® FiltrEX™ 96 and 384 Well Filter Plates

**CORNING**  
Discovering Beyond Imagination


Helen L. Krasnow, Laurence M. Vaughan<sup>2</sup>, Elizabeth M. Martin, Paul Coughlin, Linda S. Belkowski<sup>2</sup>, Debra S. Hoover  
Corning Life Sciences, Kennebunk, ME 04043, <sup>2</sup>Acton, MA 01720

## Abstract


The industrialization of genomics and proteomics demands the miniaturization of filtration devices to a more automation-compatible format. In response to this need, Corning has developed the FiltrEX 96 and 384 well filter plate platforms. These products are based on unique patented technology that allows individual filter disks to be integrally sealed in a microplate format. The objective of the current study was to demonstrate the utility of these filter plates in high throughput genomic and proteomic sample preparation methods using standard vacuum and centrifugation techniques. High throughput methods were developed for plasmid, PCR fragment and sequencing reaction clean up as well as protein size separation and desalting/concentration. The results clearly demonstrate that these products can be used to perform rapid, low cost, high quality purification of both DNA and protein and that they are compatible with standard laboratory equipment and automation.



**Figure 1. FiltrEX 96 well filter plate design**  
Proprietary nozzle design and individually integrally sealed filter disks prevents filtrate cross contamination. Rigid design and wide skirt allows for robotic handling and bar coding. **US Patent 6,391,241**




**Figure 3. Absence of cross-talk in FiltrEX 384 well filter plate**  
DNA and control buffer samples were purified in alternating wells of a FiltrEX 384 PVDF filter plate. Samples were eluted into a receiver plate and amplified by PCR. Note the absence of contaminating product in control wells.




**Figure 5. Protein Size Exclusion Using 100 K MWCO UF membrane in a 96 well plate format.**  
Protein ladder was filtered through a 100 K molecular weight cut off (MWCO) membrane and the filtrates and retentates were collected. Note exclusion of 185 kD protein.

## FiltrEX Methodology Overview



Detailed protocols available at [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or e-mail us at [genomics@corning.com](mailto:genomics@corning.com) for up to date information



**Figure 2. Plasmid and PCR samples purified using FiltrEX 96 well glass fiber plate.**  
Lanes 1, 2 and 3 contain plasmid DNA. Lanes 4, 5 and 6 contain purified PCR reactions. Lane 7 contains the unpurified PCR reaction.



**Figure 4. FiltrEX Seq UF cleanup of Big Dye® sequencing reaction.**  
ABI Version 3.1, 4 µl undiluted  
**Q20 Phred = 829**

ABI 3730: Big Dye (version 1.0)						Phred
Clone	Process	Big Dye - µl	Resusp. vol.	Q20	Classification	
FiltrEX Seq 1	10 min / 2,000 x g	2	40 µl	815	good	
FiltrEX Seq 2	10 min / 2,000 x g	2	40 µl	735	good	
FiltrEX Seq 3	10 min / 2,000 x g	2	40 µl	828	good	
FiltrEX Seq 4	10 min / 2,000 x g	2	30 µl	845	good	
FiltrEX Seq 5	10 min / 2,000 x g	2	30 µl	823	good	
FiltrEX Seq 6	10 min / 2,000 x g	2	30 µl	829	good	
Sephadex 1	2 x Wash + 5 min / 2,000 x g	2	NA	816	good	
Sephadex 2	2 x Wash + 5 min / 2,000 x g	2	NA	823	good	
Sephadex 3	2 x Wash + 5 min / 2,000 x g	2	NA	783	good	

ABI 3700: Big Dye (version 3.1)						Phred
Clone	Process	Big Dye - µl	Resusp. vol.	Q20	Classification	
FiltrEX Seq 1	10 min / 3,000 x g	1	30 µl	819	good	
FiltrEX Seq 2	10 min / 3,000 x g	2	30 µl	739	good	
FiltrEX Seq 3	10 min / 3,000 x g	4	30 µl	747	good	
Sephadex 1	2 x Wash + 5 min / 990 x g	1	NA	704	good	
Sephadex 2	2 x Wash + 5 min / 990 x g	2	NA	786	good	
Sephadex 3	2 x Wash + 5 min / 990 x g	4	NA	806	good	

**Table 1. Sequence data from FiltrEX Seq 96 vs. Sephadex G-50 cleanup.**  
Big Dye sequencing reactions were purified by FiltrEX Seq 96 (UF) or Sephadex G-50. The results are equivalent to Sephadex G-50 cleanup. Sequence was generated by Lark, Inc. (top panel) and GATC Biotech, AG (bottom panel)

## Conclusions

- FiltrEX 96 & 384 well filter plates incorporate design features that ensure well integrity and prevent sample cross contamination
- The filter plates are compatible with both vacuum and centrifugation protocols.
- The rigid skirted design meets SBS standards and ensures bar code and automation compatibility.
- FiltrEX Seq plates offer a simpler alternative to gel filtration methods for dye terminator removal.
- The availability of different media types allows the user to perform a wide variety of common genomic and proteomic sample prep methods in an easy to use, standardized format.