

Corning® Elplasia® 12K Flask

Culture Thousands of Spheroids in a Convenient Flask Format

CORNING

With the effectiveness of 3D cell culture in many areas of research, including anti-cancer drug screening and *in vitro* tumor studies, the need for better methods to produce replicate spheroids of uniform size in mass quantities has emerged.

The Corning Elplasia 12K flask addresses this need by enabling researchers to generate a high density of spheroids in a scaffold-free model.

The Corning Elplasia 12K flask contains 152 microcavities per cm² in a vessel footprint similar to that of a T-75 flask. Gravity, in conjunction with the Corning Ultra-Low Attachment (ULA) surface, and a rounded microcavity geometry enable formation of approximately 12,000 spheroids of similar shape and size.

The ULA surface, a proprietary, animal-free, covalently bonded hydrogel, promotes the formation and easy harvesting of spheroids.

The microcavity geometry allows spheroids to remain in place during medium exchange without compromising full recovery at harvest time. The flask's internal diverter feature allows for minimal disruption of spheroids during liquid handling steps.

The Corning Elplasia 12K flask is compatible with many tumor, normal, and primary cell types often used for 3D cell culture, and may be used across many applications including:

- ▶ Drug screening
- ▶ Cancer/Tumor biology
- ▶ Stem cell biology
- ▶ Cell therapy research
- ▶ 3D tissue engineering


Key Features

- ▶ Corning ULA surface
- ▶ Optically clear, gas permeable polystyrene film with low autofluorescence
- ▶ Liquid diverter feature to minimize impact of liquid flow during medium exchange steps and flask handling
- ▶ 80 cm² surface area with 152 microcavities per cm² generate approximately 12,000 spheroids of uniform size and shape per flask
- ▶ Round well-bottom geometry
- ▶ Microcavity geometry dimensions: 850 x 650 μm (top diameter x depth), with working spheroid growth dimensions of 500 x 600 μm (diameter x depth)
- ▶ 25 to 50 mL working volume

Key Benefits

- ▶ Ease of spheroid formation, culture, assessment, and harvest
- ▶ Creates uniform spheroid formation at large volumes with an easy to use “plug and play” ready protocol
- ▶ Scaffold-free cultures
- ▶ Culture spheroids for up to 30 or more days (cell line-dependent)
- ▶ Highly reproducible bulk spheroid formation across microcavities
- ▶ Common media reservoir for equivalent culture conditions for all spheroids
- ▶ Compatible with brightfield and fluorescent microscopy







The Corning® Elplasia® 12K flask was designed to enable the generation of a large volume of uniform spheroids under one culture condition with a simple “plug and play” protocol. For more details, please refer to the Corning Elplasia 12K flask Guidelines for Use (CLS-AN-713DOC).

Table 1. Cell lines evaluated using the Corning Elplasia 12K flask.


Cell Type	Source
Cancer Lines	
MCF7	Human breast adenocarcinoma
HT-29/GFP	Human colon adenocarcinoma
DU-145	Human prostate carcinoma
A549/GFP	Human lung carcinoma
Stem Cells	
hMSC	Human bone marrow mesenchymal stem cells
iPSC	Human induced pluripotent stem cells, normal
Normal Lines	
HEK-293/RFP	Human embryonic kidney
HUVEC	Human umbilical vein




The Corning Elplasia 12K flask was designed to enable the generation of large quantities of uniform spheroids. 14-day old MCF 7 (human breast adenocarcinoma) spheroids culture in a Corning Elplasia 12K flask (A). Spheroid collection post harvest in a petri dish (B). Representative image of harvested MCF 7 spheroids (C). Micrographs were taken with an EVOS® FL microscope in brightfield mode using a 2X objective.



The Corning Elplasia 12K flask substrate contains approximately 12,000 round bottom cavities that each provide a growth area of 500 x 600 μm (diameter x depth) with gas permeable bottoms.



The Corning Elplasia 12K flask is compatible with fluorescent imaging. 30-day HepG2 (human liver carcinoma) spheroids stained with Calcein AM (A), Hoechst (B), and Propidium Iodide (C). A, B, and C superimposed under Brightfield (D). Micrographs were taken with an EVOS FL microscope using a 2X objective.



The Corning Elplasia 12K flask enables long-term cultures. A-549/GFP (human lung carcinoma) spheroids cultured in a Corning Elplasia 12K flask for 30 days. Micrographs were taken with an EVOS FL microscope in brightfield and GFP modes using a 2X objective.

Ordering Information

Product may not be available in all markets.

Cat. No.	Description	Approximate Spheroids/Flask	Microcavity Size (Diameter x Depth)	Qty/Pk	Qty/Cs
4537	Corning® Elplasia® 12K flask, ULA surface, sterile	12,000	850 x 650 µm	1	5

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences
www.corning.com/lifesciences

NORTH AMERICA
t 800.492.1110
t 978.442.2200

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India
t 91 124 4604000
Japan
t 81 3-3586 1996
Korea
t 82 2-796-9500
Singapore
t 65 6572-9740
Taiwan
t 886 2-2716-0338

EUROPE
CSEurope@corning.com
France
t 0800 916 882
Germany
t 0800 101 1153
The Netherlands
t 020 655 79 28
United Kingdom
t 0800 376 8660

All Other European Countries
t +31 (0) 206 59 60 51

LATIN AMERICA
grupoLA@corning.com
Brazil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400