# Corning<sup>®</sup> Nu-Serum<sup>™</sup> and Nu-Serum IV

Growth medium supplements for reduced-serum conditions

## CORNING

#### Introduction

Reduced serum and serum-free cell culture conditions are rapidly becoming requirements for state of the art research. Corning Nu-Serum growth medium supplements provide low protein alternatives to newborn calf, fetal bovine, and other sera routinely used for cell and tissue culture.

While serum played an irreplaceable role in the early development of the ability to culture cells and tissues, it is now known that the high protein content of serum can complicate protein purification, virus production, purification and concentration, and the production and screening of monoclonal antibodies.

In comparison to standard Fetal Bovine Serum (FBS), Nu-Serum growth medium supplements are a more cost-effective method, providing the same results as a 1:1 replacement for FBS. Nu-Serum replacements have been used successfully on a large variety of human and animal cell types, many of which were previously difficult to grow.

#### **Nu-Serum and Nu-Serum IV Features**

- Full volume replacement for FBS
- Cost-effective alternative to FBS
- Facilitates protein purification, virus production, monoclonal antibody product and screening, similar to classical FBS. It increases the frequency of successful transfection of cells.
- Filtered (0.2 μm membrane)
- Supports a large variety of human and animal cell types including:
  - Human and chick embryo fibroblasts<sup>1,2</sup>
  - HeLa cells<sup>3</sup>
  - Mouse L cells<sup>4</sup>
  - HepG2 Human hepatocellular carcinoma cells⁵
  - BALB/c3T3 cells<sup>6</sup>
  - Cos cells<sup>7</sup> normal rat kidney cells<sup>8</sup>
  - Human respiratory epithelial cells9
  - Rat and chick neurons<sup>10,11</sup>
  - Osteoblasts and primary chrondrocytes<sup>12</sup>
  - Human bladder carcinoma cells13
  - Human melanomas<sup>14</sup>
  - Hybridomas derived from Sp/2 and NS-1 myelomas<sup>15</sup>





#### Potential savings using Corning<sup>®</sup> Nu-Serum<sup>™</sup> as a replacement for FBS

Medium w/10% FBS		Medium w/5% FBS		Medium w/10% Nu-Serum		Medium w/5% Nu-Serum	
500 mL medium	\$15.00	500 mL medium	\$15.00	500 mL medium	\$15.00	500 mL medium	\$15.00
55 mL FBS	\$27.50	27.5 mL FBS	\$13.50	55 mL Nu-Serum	\$20.90	27.5 mL Nu-Serum	\$10.45
Total cost	\$42.50	Total cost	\$28.50	Total cost	\$35.90	Total cost	\$25.45

Assumptions: Medium costs of \$15/500 mL bottle; premium fetal bovine serum from USDA approved countries costs of \$250/500 mL bottle. Each time a 500 mL bottle is used with a range of 5-10% supplementation, there are significant savings to the customer.

### **Ordering Information**

Corning Nu-Serum Growth Medium Supplement	25% Newborn Calf Serum		
	25% NewDorn Call Serum	100 mL	1
Corning Nu-Serum Growth Medium Supplement	25% Newborn Calf Serum	500 mL	1
Corning Nu-Serum IV Growth Medium Supplement	25% Fetal Bovine Serum	100 mL	1
Corning Nu-Serum IV Growth Medium Supplement	25% Fetal Bovine Serum	500 mL	1
	Corning Nu-Serum Growth Medium Supplement Corning Nu-Serum IV Growth Medium Supplement	Corning Nu-Serum Growth Medium Supplement25% Newborn Calf SerumCorning Nu-Serum IV Growth Medium Supplement25% Fetal Bovine Serum	Corning Nu-Serum Growth Medium Supplement25% Newborn Calf Serum500 mLCorning Nu-Serum IV Growth Medium Supplement25% Fetal Bovine Serum100 mL

#### References

- 1. Davis, M.H., Arch Biochem. Biophys., 251:498 (1987)
- 2. Goodman,S.L., et al., Cell, 41:1029 (1985)
- 3. Medh, R.D., et al., Blood, 80:981 (1992)
- 4. Chu, Y.-W., et al., Proc. Nat. Acad. Sci. USA, 20:4261 (1993)
- 5. Lucore, D.C., et al., J. Biol. Chem, 263:15845 (1988)
- 6. Guity, M., et al., J. Cell Biol., 116:1055 (1992)
- 7. Attisano, L., et al., Cell, 68:97 (1992)
- 8. Cluett, E.B., et al., J. Cell Biol., 120:15 (1993)
- 9. Devalia, J.L., et al., Resp. Med., 84:303 (1990)
- 10. Rage, F., et al., Endrocroniol., 130:1056 (1992)
- 11. Hory-Lee, F., et al., Proc. Nat. Acad. Sci USA, 90:2613 (1993)
- 12. Elfort, P.R., et al., Endrocrinol, 127:1635 (1990)
- 13. Guity, M., et al., J. Cell Biol., 111:2765 (1990)
- 14. Matsuo, 0., et al., Acta Haemotol, Japan, 47:1049 (1984)
- 15. Sjogren-Jansson, E., and Jeansson, S., J. Immunol. Meth., 84:359 (1985)

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

#### Mediatech, Inc.

A Corning Subsidiary 9345 Discovery Boulevard Manassas, VA 20109 t 800.235.5476 t 703.471.5955 f 703.467.9851

www.corning.com/lifesciences/ media For additional product or technical information, visit **www.corning.com/lifesciences** or call 1.800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

GOSSELIN

**PYREX** 

AXYGEN

### CORNING |

For a listing of Corning trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.

FALCON