

# Corning® hybrigro SF™ Medium

## Animal-Free, Serum-Free Defined Medium

### with Corning glutagro™ Supplement

CORNING  
cellgro®




### Introducing Corning hybrigro SF Medium

hybrigro SF Medium is a complete, animal component-free, defined medium specifically developed for serum-free growth and antibody production with a variety of hybridoma cell lines.

### Features

- ▶ Improved growth and antibody production over serum supplemented media
- ▶ Formulated with glutagro supplement to provide a stable, ready-to-use formulation; no glutamine supplementation required
- ▶ Contains Polaxamer 188™ to provide shear force protection during agitated suspension culture
- ▶ Formulated without phenol red.
- ▶ Contains low concentrations of highly purified non-animal origin recombinant proteins added for enhanced cell growth and productivity
- ▶ Defined: there are no undefined hydrolysates, ultra filtered lysates, or peptides in the formulation



### Results

Comparison of hybrigro SF Medium versus standard serum-supplemented medium using AE1 (ATCC HB-72™) hybridoma cell line. Growth is superior with hybrigro SF medium, with enhanced viability.

Note: hybrigro SF Medium will not grow cholesterol-auxotrophic myeloma-derived hybridomas (e.g., NS0, NS1, and P3X63Ag8.653) without further supplementation. Supplementation with a solubilized preparation of cholesterol is required for these cell line derivatives.

## Corning® hybrigro™ Medium Ordering Information

Corning cellgro Cat. No.	Description	Unit Size	Qty/Pk
40-215-CV	hybrigro SF Medium with glutagro supplement	500 mL	6

For additional product or technical information, please visit [www.cellgro.com](http://www.cellgro.com) or call 1.800.235.5476.

**Mediatech, Inc.**  
*A Corning Subsidiary*  
9345 Discovery Boulevard

Manassas, VA 20109  
t 800.235.5476  
t 703.471.5955  
f 703.467.9851

[www.cellgro.com](http://www.cellgro.com)

*The Corning Family of Brands*

CORNING costar cellgro PYREX AXYGEN GOSSELIN™

For a listing of Corning trademarks, visit us at [www.cellgro.com/about-us/trademarks.html](http://www.cellgro.com/about-us/trademarks.html)  
All other trademarks in this document are the property of their respective owners.  
Corning Incorporated, One Riverfront Plaza, Corning, NY 14831-0001