

# Corning® Synthege<sup>™</sup> 3D hiPSC Suspension Matrix Kit

CORNING

## Guidelines for Use

The Corning Synthege hiPSC Suspension Matrix Kit is a powerful bio-tool for large-scale manufacturing of physiological 3D human induced Pluripotent Stem Cells (hiPSCs) spheroids in a lab setting. Cell encapsulation and spheroid isolation is an easy and straightforward process. For example: one 6-well plate seeded with hiPSC at a density of  $5 \times 10^4$  cell/mL cultured for 5 days using the Corning Synthege hiPSC Suspension Matrix Kit, can yield approx. 50 million total cells (approx. 600,000 to 700,000 spheroids with diameter ranging from 30  $\mu\text{m}$  to 50  $\mu\text{m}$ ). Additionally, daily media exchanges are not required to maintain hiPSC as in more traditional planar culture. After retrieving, spheroids can be used directly for various downstream applications (e.g., drug screening, bioprinting for tissue engineering, and somatic cells differentiation). The Corning Synthege hiPSC Suspension Matrix Kit consists of a vial of Corning Synthege hiPSC proprietary peptide solution for 3D suspension culture and a vial of Corning Synthege X-Link solution. The Synthege hiPSC Suspension Matrix nanofibrils are formulated into a basal cell culture medium containing Corning Synthege 3D hiPSC Grow Mix. This allows for a neutral pH that forms a 3D microenvironment suitable for spheroid growth. With Corning Synthege hiPSC Suspension Matrix Kit, cells no longer suffer acidic or chill conditions since all operating procedures can be completed at room temperature or 37°C and at neutral pH. Additionally, there is no need to stir the cells as typically required in suspension culture, so cells are not exposed to shear force. There is also flexibility on the volume of Synthege hiPSC Suspension Matrix Kit used per vessel which allows users to scale-up as needed utilizing any standard cell culture vessel.

**NOTE:** In addition to the Corning Synthege hiPSC Suspension Matrix Kit, Corning Synthege 3D hiPSC Grow Mix (Cat. No. 354792) needs to be purchased separately.

### 3D Suspension hiPSC Culture Protocol

Bring the Synthege hiPSC Suspension Matrix Kit (Synthege hiPSC Matrix Peptide solution and Synthege X-Link solution) and culture medium to room temperature.

#### Step 1. Complete culture medium stock solution preparation

1. Add 300  $\mu\text{L}$  sterile DPBS (without  $\text{Mg}^{2+}$  and  $\text{Ca}^{2+}$ ) into one (1 mg) vial of lyophilized Synthege hiPSC Grow Mix. Pipette gently to obtain a homogeneous solution.
  - The reconstituted Synthege hiPSC Grow Mix solution is stable at  $-20^\circ\text{C}$  for up to 6 months, therefore, aliquots should be prepared, with working volumes, to avoid repeated freeze-thaw cycles.
2. Add an aliquot of the reconstituted Synthege hiPSC Grow Mix solution into pre-warmed cell culture medium (i.e., mTeSR<sup>™</sup>1 complete medium) at a ratio of 1:1000 v/v (Synthege hiPSC Grow Mix solution: mTeSR1 complete medium) to prepare hiPSC complete medium.
  - All medium used for cell culture in this protocol needs Synthege hiPSC Grow Mix and should be used within 2 weeks of addition.

Synthege hiPSC Grow Mix contains ROCK Inhibitor and does not require additional ROCK supplementation.

#### Step 2. Prepare Synthege hiPSC Suspension Matrix Kit

**NOTE:** Corning Synthege hiPSC Suspension Matrix Kit can be used at volumes between 0.3 and 0.85 mL/cm<sup>2</sup>. This allows the volume to be adjusted depending on how many cells are needed. See table below for an example of a culture using 0.3 or 0.85 mL/cm<sup>2</sup>.

Corning Synthege Thickness (mL/cm <sup>2</sup> )	Volume of Corning Synthege hiPSC Suspension Matrix Kit per T-75 Flask	Cells Required to Seed (assuming 50,000 cells/mL)	Cells at Harvest (assuming 3 population doublings)
0.3	22.5 mL	$1.125 \times 10^6$ cells	$9 \times 10^6$ cells
0.85	63.75 mL	$3.15 \times 10^6$ cells	$2.52 \times 10^7$ cells

1. The table below can be used to calculate number of cells, media, X-Link solution and Synthege matrix (peptide) solution needed based on the listed % of total volume determined above.

- We recommend starting with a cell density of between 40-100,000 cells/mL.

Media/Cells	Corning X-Link Solution	Corning Synthege hiPSC Suspension Matrix (Peptide)
84%	0.5%	15.5%

**Example:** A single well of a 6-well plate at 0.85 mL/cm<sup>2</sup>

Total volume of Synthegel hiPSC Suspension Matrix Kit at 0.85 mL/cm<sup>2</sup>

- 0.85 mL/cm<sup>2</sup> x 9.5 cm<sup>2</sup> (surface area) = 8.075 mL

Volume of each component of Synthegel hiPSC Suspension Matrix

- Media/cells: 84% of 8.075 mL = 6.80 mL
- Cells (i.e., desired cell density of 50,000 cells/mL): 8.075 mL x 50,000 cells/mL = 4 x 10<sup>5</sup> cells in 6.80 mL media
- X-Link Solution: 0.5% of 8.075 mL = 40 µL
- Synthegel Peptide: 15.5% of 8.075 mL = 1.25 mL

2. Prepare 6.80 mL cell suspension Mixture A (Figure 1) containing 4 x 10<sup>5</sup> cells using hiPSC complete medium.

3. Cell Encapsulation: Add 1.25 mL Synthegel peptide solution to Mixture A. Mix gently with a pipet to avoid making bubbles.

4. Add 40 µL Synthegel X-Link solution to Mixture B resulting in the formation of Mixture C with a final cell density 5 x 10<sup>4</sup> cell/mL. Mix gently with a pipet to avoid making bubbles.

5. Transfer 8.09 mL of Mixture C to a single well of a 6-well-plate.


### Step 3: Feeding strategy

**NOTE:** Unlike iPSCs grown on two-dimensional formats, iPSC grown in the hiPSC Suspension Matrix Kit may not require daily media exchanges.

1. Incubate the plate at 37°C (5% CO<sub>2</sub>) for 4 to 6 days. Additional complete media containing Grow Mix can be added on Day 4 or 5 of culture by gently pipetting or swirling.

- Do not add more than 30% of the initial Synthegel hiPSC Suspension Matrix Kit volume to culture as too much media can dilute the suspension matrix and result in undesired cell sedimentation and or attachment.

**Example:** 24% of 8.075 mL = 1.9 mL total of media. Feed media can be added on a single day (Day 3, 4 or 5) or divided across Days 3 to 5.



**Figure 1.** Corning Synthegel hiPSC Suspension Matrix Kit cell culture (example for 6-well plate).

### Step 4. Harvest of hiPSC spheroids

**Example:** A single well of a 6-well plate

1. Mechanically disrupt the 3D suspension culture by pipetting 5 to 6 times in the well. Transfer to a 15 mL centrifuge tube.
2. Use 2 mL Dulbecco's Phosphate-Buffered Saline (DPBS, without Mg<sup>2+</sup> and Ca<sup>2+</sup>) to rinse the well and transfer to the centrifuge tube.
3. Centrifuge at 500 g for 5 minutes by using a centrifuge equipped with a swing bucket rotor. Discard the supernatant, and collect the pellet containing cell spheroids.


## Step 5. Cell dissociation

**NOTE:** For dissociation of spheroids, Accutase® cell detachment solution (Corning) or 10X TrypLE™ is recommended.

1. Add 6 to 7 mL Accutase cell detachment solution to the cell pellet, and pipet several times to resuspend pellet.
2. Incubate at room temperature for 15 to 20 min. Visualization under the microscope can help ensure spheroid dissociation.
3. Pipette gently to help breakup the spheroids to desired size.
4. Process cells as needed.

To start a culture from freeze-thaw, a higher starting cell seeding density is recommended.

- Thaw of 3D cultured cells:  $1-2 \times 10^5$  cell/mL.
- Thaw of 2D cultured cells:  $2-4 \times 10^5$  cell/mL



**Figure 2.** Morphologies of hiPSC spheroids after 5 days of culture in Corning SyntheGel hiPSC Suspension Matrix Kit with mTeSR Plus medium. Initial seeding density was  $4.5 \times 10^4$  cell/mL.

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use or general laboratory use only.\* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. \*For a listing of US medical devices, regulatory classifications or specific information on claims, visit [www.corning.com/resources](http://www.corning.com/resources).

*Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.*

# CORNING

Corning Incorporated  
Life Sciences  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

**NORTH AMERICA**  
t 800.492.1110  
t 978.442.2200

**ASIA/PACIFIC**  
**Australia/New Zealand**  
t 61 427286832  
**Chinese Mainland**  
t 86 21 3338 4338

**India**  
t 91 124 4604000  
**Japan**  
t 81 3-3586 1996  
**Korea**  
t 82 2-796-9500  
**Singapore**  
t 65 6572-9740  
**Taiwan**  
t 886 2-2716-0338

**EUROPE**  
CSEurope@corning.com  
**France**  
t 0800 916 882  
**Germany**  
t 0800 101 1153  
**The Netherlands**  
t 020 655 79 28  
**United Kingdom**  
t 0800 376 8660

**All Other European Countries**  
t +31 (0) 206 59 60 51

**LATIN AMERICA**  
grupoLA@corning.com  
**Brazil**  
t 55 (11) 3089-7400  
**Mexico**  
t (52-81) 8158-8400