

# The Corning® 5L Erlenmeyer Flask Provides an Ideal Choice for Efficient Scale-up of Suspension Cells

CORNING

## SnAPPShots

A brief technical report from the Corning Applications Group

Hilary Sherman, B.S.  
Corning Incorporated, Life Sciences  
Kennebunk, Maine 04043

### Introduction


As the bioprocess industry grows, there is a corresponding growth in demand for single-use containers that can efficiently culture large quantities of cells. To aid in this effort, Corning has developed a 5L Erlenmeyer flask with the same footprint as a traditional 3L Erlenmeyer flask. The 5L shape has been optimized for increased gas exchange compared to the more traditional Erlenmeyer flask designs. Due to gas exchange limitations, most traditionally shaped Erlenmeyer flasks can only accommodate approximately one third of the stated volume of the vessel. Corning's new vessel design allows for the culture volume to be increased to greater than one half of the stated volume of the vessel without compromising cell viability or density. The increased volume ratio results in a greater number of cells (higher cell yield) cultured in less space. In this study, we demonstrate that the Corning 5L Erlenmeyer flask exhibits enhanced performance compared to a competitor's 5L Erlenmeyer flask when culturing a suspension cell line commonly used for baculovirus production.

### Materials and Methods


Sf9 cells (Life Technologies Cat. No. 12659-017) were thawed and cultured in Sf900™ II SFM (Life Technologies Cat. No. 10902104) at a density of  $1 \times 10^6$  to  $8 \times 10^6$  cells/mL in a Gallenkamp floor model orbital shaker at 27°C. After cell expansion, Sf9 cells were seeded into either a Corning (Cat. No. 431684) or competitor 5L Erlenmeyer flask at a density of  $1 \times 10^6$  cells/mL. Cells in 2.5L of medium were rotated at 90 rpm while cells in 3.5L of medium were rotated at 100 rpm. Cells were counted daily on the Nova BioProfile® Flex analyzer for 4 days. These studies were repeated one additional time for a total of 4 vessels per condition for 2.5L studies and 2 vessels per condition for 3.5L studies.

### Results

Throughout the 4 day studies, the viability of Sf9 cultures were maintained above 98%, irrespective of the vessel or rpm that was utilized (Figures 1 and 2). After 4 days of growth at 2.5L of medium, the competitor vessel exhibited statistically lower cell



**Figure 1.** Sf9 viability was maintained above 98% regardless of which vessel was used (right axis). After 3 days of growth in 2.5L, the competitor vessel exhibited lower cell densities compared to the Corning 5L Erlenmeyer flask (left axis). Two way ANOVA \* =  $p < 0.05$ , \*\*\* $p < 0.001$ ,  $n = 4$ .



**Figure 2.** Sf9 viability was maintained above 98% regardless of which vessel was used (right axis). After 4 days of growth in 3.5L, the competitor vessel exhibited statistically lower cell densities compared to the Corning 5L Erlenmeyer flask (left axis). Two way ANOVA \*\*\* $p < 0.001$ ,  $n = 2$ .

densities compared to the Corning® 5L Erlenmeyer flask (Figure 1). When the volume was increased to 3.5L of medium, the Corning vessel performed similarly to that at 2.5L of volume and statistically out-performed the competitor 5L flask at 3.5L of culture medium (Figure 2).


### Conclusions



- ▶ The Corning 5L Erlenmeyer flask can be used to generate greater cell yields when compared to a traditional 3L Erlenmeyer flask utilizing the same footprint.
- ▶ Sf9 cell viability was equivalent when cells were cultured in the Corning or competitor 5L Erlenmeyer flasks.
- ▶ After 3 days of growth, Corning 5L Erlenmeyer flasks exhibited statistically higher cell yields compared to the competitor equivalent when 2.5L of medium was used.
- ▶ Similar cell densities were obtained in the Corning 5L Erlenmeyer flasks when cells were cultured in 2.5L or 3.5L of culture medium.
- ▶ Corning 5L Erlenmeyer flasks exhibited statistically higher cell yields compared to the competitor equivalent when 3.5L of medium was used.

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.



At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, please call 800.492.1110 or visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences). Customers outside the United States, call +1.978.442.2200 or contact your local Corning sales office listed below.

#### Corning Incorporated Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

#### Worldwide Support Offices

**ASIA/PACIFIC**  
**Australia/New Zealand**  
t 0402-794-347

**China**  
t 86 21 2215 2888  
f 86 21 6215 2988

**India**  
t 91 124 4604000  
f 91 124 4604099

**Japan**  
t 81 3-3586 1996  
f 81 3-3586 1291

**Korea**  
t 82 2-796-9500  
f 82 2-796-9300

**Singapore**  
t 65 6733-6511  
f 65 6861-2913

**Taiwan**  
t 886 2-2716-0338  
f 886 2-2516-7500

#### EUROPE

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

#### All Other European Countries

t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

#### LATIN AMERICA

**Brasil**  
t (55-11) 3089-7419  
f (55-11) 3167-0700

**Mexico**  
t (52-81) 8158-8400  
f (52-81) 8313-8589

**CORNING** | **FALCON** | **AXYGEN** | **GOSELIN** | **PYREX**