

# P251 Efficient Solute Carrier (SLC) Transporter-Mediated Drug-Drug Interaction Testing in a 96-well Plate Format

Lisa G. Fox, Molly McDonnell, Niele Olivier, Inese Smukste, **Ritu Singh**, and David M. Stresser


Corning® Gentest<sup>SM</sup> Contract Research Services, Corning Life Sciences, Corning Incorporated  
6 Henshaw Street, Woburn MA 01801, USA

CORNING

## Abstract

Membrane transporters can be major determinants of pharmacokinetic properties of drugs and their metabolites, and can mediate drug-drug interactions (DDIs). We have established an *in-vitro* testing platform for systematic analysis of solute carrier (SLC) transporter interactions, for both screening and definitive DDI studies recommended by regulatory agencies. Corning® TransportoCells™ are HEK-293 cells that transiently over-express a single SLC transporter protein and provide a common cell-based model. The transporters examined were OATP1B1, OATP1B3, OAT1, OAT3, OCT1, OCT2, MATE1 and MATE2-K, and assays were established using prototypical and/or clinically relevant substrates and inhibitors. TransportoCells were thawed, seeded in 24- and 96-well plates coated with poly-D-lysine, and cultured for 24h. Before assay, cells were washed with pre-warmed HBSS buffer (with Ca<sup>2+</sup> and Mg<sup>2+</sup>) and pre-incubated in HBSS buffer for 10 min at 37° C. MATE1 and MATE2-K cells were subsequently pre-treated with 40 mM NH<sub>4</sub>Cl in HBSS for 20 min. The uptake assays were initiated by addition of radiolabeled substrate and non-radiolabeled inhibitor, if applicable, and incubated at 37° C. Reactions were stopped by removing the dosing solutions and washing cells with cold HBSS. The cells were lysed with M-Per Mammalian Protein extraction reagent, and the uptake activity was quantified using LSC normalized for protein concentration in each sample. K<sub>m</sub> and IC<sub>50</sub> values for OATP1B1, OATP1B3, OAT1, OAT3, OCT1, OCT2, MATE1 and MATE2-K are shown in the table below. We found robust performance in a 96-well plate format with only minimum medium volumes (to conserve consumption of certain costly radiolabeled substrates). Efficiency was further improved by eliminating use of sodium butyrate pretreatment for OATP1B1 and OATP1B3 for prototypical substrates. We have established a robust *in-vitro* testing system based on a common cell model of HEK-293 cells transiently over-expressing a single SLC transporter, for conducting screening and definitive transporter-mediated DDI studies recommended by regulatory agencies.

Figure 1. Expression of clinically important hepatic and renal transporters



## Materials and Methods

**Materials:** Corning® TransportoCells™ products: OATP1B1\*1a (Cat. No. 354859), OATP1B3 (Cat. No. 354851), OAT1 (Cat. No. 354857), OAT3 (Cat. No. 354858), OCT1 (Cat. No. 354852), OCT2 (Cat. No. 354853), MATE1 (Cat. No. 354855), MATE2-K (Cat. No. 354856), and empty vector (control) cells (Cat. No. 354854) and cell culture reagents were obtained from Corning Life Sciences. Radiolabeled and non-radiolabeled substrates and inhibitors: estradiol-17 beta-D-glucuronide (E<sub>2</sub>-17-β-Gluc), cholecystokinin octapeptide (CCK-8), estrone-3-sulfate (E3S), p-amino hippuric acid (PAH), metformin, rifampicin, probenecid, and cimetidine were obtained from Corning, PerkinElmer, American Radiolabeled Chemicals, or Moravak.

**Inhibition Assay:** Corning® TransportoCells™ products transiently over-expressing a single SLC transporter and control cells transfected with empty vector were thawed, seeded in 24-well or 96-well plates coated with poly-D-lysine, and cultured for 24 h. Before assay, cells were washed with pre-warmed HBSS buffer (with Ca<sup>2+</sup> and Mg<sup>2+</sup>) and pre-incubated in HBSS buffer for 10 min at 37° C. MATE1 and MATE2-K were subsequently pre-treated with 40 mM NH<sub>4</sub>Cl in HBSS for 20 min. The uptake assays were initiated by addition of radiolabeled substrate and non-radiolabeled inhibitors, if applicable, and incubated at 37° C for 2 min (MATE1/2K), 5 min (OATP1B1\*1a, OATP1B3, OAT1/3), or 10 min (OCT1/2). Reactions were stopped by removing the dosing solutions and washing the cells with cold HBSS. The cells were lysed with M-Per Mammalian Protein extraction reagent, and the uptake activity was quantified using liquid scintillation counting normalized for protein concentration in each sample.

**Data Analysis:** Kinetic parameters were determined by non-linear regression using XLfit (IDBS). For each substrate concentration, the initial uptake rate was calculated by subtracting the initial rate determined in HEK-293 cells expressing empty vector from those obtained in HEK-293 cells over-expressing a SLC transporter. For the inhibition assays, IC<sub>50</sub> values were determined by using Sigmoidal Hill four-parameter equation.

## Results and Discussion

Initially, we established the Corning TransportoCells assay for contract research services using the 24-well format, obtaining K<sub>m</sub>, V<sub>max</sub>, and IC<sub>50</sub> results with prototypical substrates, which are consistent with reported values in the literature. Assay conditions were then modified to incorporate operational efficiencies while retaining robust performance:

- The assay was established in the 96-well plate format, allowing volume reduction of TransportoCells product, medium, and certain costly radiolabeled substrates.
- The savings in reagent volumes using 96-well format vs. 24-well, were further increased by reducing the well volumes 33% from 120 µL to 80 µL.

- Sodium butyrate can be added to culture medium to enhance transporter activity. We found that omission of sodium butyrate from the medium (OATP1B1 and OATP1B3 assays with prototypical substrates) resulted with adequate activity of the substrates. Moreover, we observed an enhancement of cell attachment when sodium butyrate was not added to the medium. Therefore, supplementing the medium with sodium butyrate is recommended when additional activity is required, however, the step can be omitted when signal is adequate, thereby reducing the steps in the assay.

The K<sub>m</sub>, V<sub>max</sub>, and IC<sub>50</sub> values for OATP1B1\*1a, OATP1B3, OAT1, OAT3, OCT1, OCT2, MATE1, and MATE2-K generated in the 96-well format were in good agreement with our 24-well results. Representative results of our evaluation are summarized in Table 1, Table 2, and Table 3.

Results from a representative Corning Gentest<sup>SM</sup> Contract Research study performed in the 96-well Corning TransportoCells assay are shown in Table 5.

Table 1. Comparison of uptake (K<sub>m</sub>) results for prototypical substrates of SLC transporters in 96-well vs. 24-well Corning TransportoCell assays

| Transporter | Substrate                 | Plate Format | K <sub>m</sub> (µM) | V <sub>max</sub> (pmol/mg/min) | Literature K <sub>m</sub> Value (µM) | Literature Reference |
|-------------|---------------------------|--------------|---------------------|--------------------------------|--------------------------------------|----------------------|
| OATP1B1*1a  | E <sub>2</sub> -17-β-Gluc | 24           | 9.8                 | 337                            | 6.3                                  | 4                    |
|             |                           | 96           | 8.5                 | 424                            |                                      |                      |
| OATP1B3     | CCK-8                     | 24           | 39                  | 4998                           | 17                                   | 5                    |
|             |                           | 96           | 52                  | 4649                           |                                      |                      |
| OAT3        | E3S                       | 24           | 3.7                 | 276                            | 6.3                                  | 6                    |
|             |                           | 96           | 8.1                 | 618                            |                                      |                      |

Table 2. Uptake (K<sub>m</sub>) data for prototypical substrates of SLC transporters in 96-well Corning TransportoCell assays

| Transporter | Substrate                 | K <sub>m</sub> (µM) |         | V <sub>max</sub> (pmol/mg/min) |         | Literature K <sub>m</sub> | Literature Reference |
|-------------|---------------------------|---------------------|---------|--------------------------------|---------|---------------------------|----------------------|
|             |                           | Assay 1             | Assay 2 | Assay 1                        | Assay 2 |                           |                      |
| OATP1B1*1a  | E <sub>2</sub> -17-β-Gluc | 6.0                 | 11      | 337                            | 511     | 6.3                       | 4                    |
| OATP1B3     | CCK-8                     | 37                  | 67      | 3277                           | 6021    | 17                        | 5                    |
| OAT1        | PAH                       | 99                  | 118     | 3941                           | 6305    | 28                        | 6                    |
| OAT3        | E3S                       | 9.1                 | 7.1     | 863                            | 373     | 6.3                       | 6                    |
| OCT1        | Metformin                 | 8165                | 5679    | 22899                          | 17287   | 5450                      | 7                    |
| OCT2        | Metformin                 | 3766                | 6618    | 52545                          | 31703   | 3356                      | 8                    |
| MATE1       | Metformin                 | 541                 | 333     | 34281                          | 22384   | 227                       | 9                    |
| MATE2-K     | Metformin                 | 2408                | 2476    | 48932                          | 29134   | 1980                      | 10                   |

Table 3. Inhibition of uptake (IC<sub>50</sub>) data for prototypical substrates of SLC transporters in 96-well Corning TransportoCell Assays

| Transporter | Substrate                 | Inhibitor  | IC <sub>50</sub> (µM) |         | Literature IC <sub>50</sub> | Literature Reference |
|-------------|---------------------------|------------|-----------------------|---------|-----------------------------|----------------------|
|             |                           |            | Assay 1               | Assay 2 |                             |                      |
| OATP1B1*1a  | E <sub>2</sub> -17-β-Gluc | Rifampicin | 1.4                   | 1.9     | 1.5                         | 11                   |
| OATP1B3     | CCK-8                     | Rifampicin | 5.0                   | 4.6     | no data                     | no data              |
| OAT1        | PAH                       | Probenecid | 6.5                   | 5.7     | 6.5                         | 12                   |
| OAT3        | E3S                       | Probenecid | 7.8                   | 7.7     | 9.0 (K)                     | 13                   |
| OCT1        | Metformin                 | Cimetidine | 109                   | 251     | 166 (K)                     | 14                   |
| OCT2        | Metformin                 | Cimetidine | 446                   | 642     | 373                         | 15                   |
| MATE1       | Metformin                 | Cimetidine | 2.2                   | 2.0     | 3.8 (K)                     | 16                   |
| MATE2-K     | Metformin                 | Cimetidine | 6.7                   | 7.4     | 6.9                         | 16                   |

In addition, we conducted further characterization of OATP1B1\*1a- and OATP1B3-expressing Corning TransportoCells products with clinically relevant substrates, statin drugs pravastatin and rosuvastatin. The results are shown in Table 4 and representative K<sub>m</sub>/IC<sub>50</sub> curves are provided in Figure 2.

Table 4. Kinetic data for statins in 96-well Corning TransportoCell assays

| Transporter | Substrate    | K <sub>m</sub> (µM) | V <sub>max</sub> (pmol/mg/min) | Inhibitor  | IC <sub>50</sub> value (µM) |
|-------------|--------------|---------------------|--------------------------------|------------|-----------------------------|
| OATP1B1*1a  | Pravastatin  | 58                  | 238                            | Rifampicin | 5.2                         |
| OATP1B1*1a  | Rosuvastatin | 11                  | 76                             | Rifampicin | 2.0                         |
| OATP1B3     | Pravastatin  | 70                  | 107                            | Rifampicin | 0.74                        |
| OATP1B3     | Rosuvastatin | 9.7                 | 140                            | Rifampicin | 14                          |

Figure 2. K<sub>m</sub> and IC<sub>50</sub> curves for statins in 96-well Corning TransportoCells assays




Table 5. Results from a representative Corning® Gentest<sup>SM</sup> Contract Research Services study (Transporter DDI) conducted in Corning TransportoCells

| Transporter | Probe Substrate                          | Inhibitor      | Inhibitor Conc. (µM) | % of Remaining Activity | Uptake Ratio |
|-------------|--|----------------|----------------------|-------------------------|--------------|
| OATP1B1     | 1 µM E <sub>2</sub> -17-β-Gluc           | Test Article X | 0                    | 100%                    | 10           |
|             |  | Rifampicin     | 300                  | 64%                     | 7.6          |
|             |  | Rifampicin     | 100                  | 0.26%                   | 1.0          |
| OATP1B3     | 2 µM cholecystokinin octapeptide (CCK-8) | Test Article X | 0                    | 100%                    | 12           |
|             |  | Rifampicin     | 300                  | 76%                     | 13           |
|             |  | Rifampicin     | 100                  | 7.8%                    | 2.5          |
| OAT1        | 15 µM PAH                                | Test Article X | 0                    | 100%                    | 22           |
|             |  | Probenecid     | 300                  | 109%                    | 20           |
|             |  | Probenecid     | 100                  | 16%                     | 3.8          |
| OAT3        | 4 µM E3S                                 | Test Article X | 0                    | 100%                    | 9.5          |
|             |  | Probenecid     | 300                  | 23%                     | 3.7          |
|             |  | Probenecid     | 100                  | 7.0%                    | 1.6          |
| OCT1        | 10 µM Metformin                          | Test Article X | 0                    | 100%                    | 11           |
|             |  | Cimetidine     | 300                  | 33%                     | 5.8          |
|             |  | Cimetidine     | 1000                 | 23%                     | 3.5          |
| OCT2        | 10 µM Metformin                          | Test Article X | 0                    | 100%                    | 20           |
|             |  | Cimetidine     | 300                  | 45%                     | 11           |
|             |  | Cimetidine     | 1000                 | 16%                     | 4.8          |

A

Data generated in an SLC Transporter DDI study for a pharmaceutical Sponsor.

The inhibition screening assay measuring transporter-mediated uptake activity of prototypical substrates with six SLC transporters, identified inhibitory potential of Test Article X with three of the transporters (A).

The concentration-dependent inhibitory effect was further characterized in the subsequent IC<sub>50</sub> assay (B).

B

| Transporter | Probe Substrate | Inhibitor      | Inhibitor Conc. (µM) | % of Remaining Activity |             |             | Mean % of Remaining Activity | IC <sub>50</sub> (µM) |
|-------------|-----------------|----------------|----------------------|-------------------------|-------------|-------------|------------------------------|-----------------------|
|             |                 |                |                      | Replicate 1             | Replicate 2 | Replicate 3 |                              |                       |
| OAT3        | 4 µM E3S        | Test Article X | 0                    | 100                     | 100         | 100         | 100                          | 42                    |
|             |                 |                | 0.41                 | 102                     | 96          | 114         | 103                          |                       |
|             |                 |                | 1.2                  | 93                      | 91          | 109         | 97                           |                       |
|             |                 |                | 3.7                  | 84                      | 85          | 101         | 89                           |                       |
|             |                 |                | 11                   | 73                      | 69          | 83          | 75                           |                       |
|             |                 |                | 33                   | 60                      | 53          | 64          | 58                           |                       |
|             |                 |                | 100                  | 29                      | 27          | 33          | 29                           |                       |
|             |                 |                | 300                  | 15                      | 13          | 11          | 13                           |                       |
|             |                 |                | 1000                 | 6.3                     | 5.7         | 8.0         | 6.6                          |                       |
|             |                 |                | 1000                 | 6.3                     | 5.7         | 8.0         | 6.6                          |                       |
| OCT1        | 10 µM Metformin | Test Article X | 0                    | 100                     | 100         | 100         | 100                          | 178                   |
|             |                 |                | 0.41                 | 113                     | 68          | 101         | 94                           |                       |
|             |                 |                | 1.2                  | 105                     | 77          | 91          | 91                           |                       |
|             |                 |                | 3.7                  | 112                     | 82          | 95          | 96                           |                       |
|             |                 |                | 11                   | 77                      | 69          | 97          | 81                           |                       |
|             |                 |                | 33                   | 93                      | 82          | 97          | 90                           |                       |
|             |                 |                | 100                  | 60                      | 74          | 68          | 67                           |                       |
|             |                 |                | 300                  | 32                      | 38          | 30          | 34                           |                       |
|             |                 |                | 1000                 | 22                      | 21          | 23          | 22                           |                       |
|             |                 |                | 1000                 | 22                      | 21          | 23          | 22                           |                       |
| OCT2        | 10 µM Metformin | Test Article X | 0                    | 100                     | 100         | 100         | 100                          | >300                  |
|             |                 |                | 0.41                 | 91                      | 92          | 87          | 90                           |                       |
|             |                 |                | 1.2                  | 96                      | 83          | 74          | 85                           |                       |
|             |                 |                | 3.7                  | 88                      | 57          | 84          | 77                           |                       |
|             |                 |                | 11                   | 98                      | 94          | 95          | 96                           |                       |
|             |                 |                | 33                   | 98                      | 94          | 101         | 98                           |                       |
|             |                 |                | 100                  | 87                      | 98          | 102         | 95                           |                       |
|             |                 |                | 300                  | 46                      | 56          | 62          | 54                           |                       |
|             |                 |                | 1000                 | 14                      | 16          | 18          | 16                           |                       |
|             |                 |                | 1000                 | 14                      | 16          | 18          | 16                           |                       |

## Conclusions

We have established a robust *in-vitro*, 96-well testing system based on a common cell model (Corning TransportoCells assay – HEK-293 cells transiently over-expressing a single SLC transporter) for conducting screening and definitive DDI studies required by regulatory agencies.

Operational efficiencies and cost savings were gained by adapting the methods to the 96-well format, which included omission of a medium supplement and reduction of expensive reagent volumes.

The method has been successfully used in the Corning Gentest<sup>SM</sup> Contract Research Services organization for both screening and DDI studies contracted by pharmaceutical Sponsors.

## References

- Glacomini, KM, et al., International Transporter Consortium (2010) Nat Rev Drug Discov 9:215-236.
- US FDA Draft Guidance for Industry: Drug Interaction Studies (2012).
- EMA Guidance on Investigation of Drug Interaction (Effective 2013).
- Sharma P, et al., Xenobiotica 2010 Jan;40(1):24-37.
- Poirier A, et al., J Pharmacol Exp Ther 2009 Dec;326(6):585-611.
- Ueo H, et al., Biochem Pharmacol. 2005 Oct 1;70(7):1104-13.
- Ahlin G, et al., Pharmacogenomics J. 2011 Dec;11(6):400-11.
- Zolk O, et al., Naunyn-Schmiedeberg Arch Pharmacol. 2009 Apr;379(4):337-48.
- Chen Y, et al., Pharmacogenomics J. 2009 Apr;9(2):127-36.
- Tanihara Y, et al., Biochem Pharmacol. 2007 Jul 15;74(2):359-71.
- Gu C, et al., Eur J Pharmacol. 2008 Apr 14;584(1):57-65.
- Ho ES, et al., Am Soc Nephrol. 2000 Mar;11(3):383-93.
- Takeda M, et al., Eur J Pharmacol. 2001 May 11;419(2-3):113-20.
- hang, L, et al., J Pharmacol Exp Ther. 1998 286:354.
- Okuda M, et al., Biochim Biophys Acta. 1999 Mar 4;1417(2):224-31.
- Itto S, et al., J Pharmacol Exp Ther. 2012 Feb;340(2):393-403