Corning® BioCoat™ Culture Slides and Falcon® Culture Slides

CORNING

Guidelines for Use

ENGLISH

Intended Use

Corning BioCoat and Falcon culture slides are used for cell and tissue culture. A plastic chamber on a glass slide permits culture in defined compartments. After culture is completed, samples can be treated with immunological labeling, or staining reagents directly in the chamber. The chamber can then be removed for microscopy. The use of culture slides eliminates the need to transfer cells or tissue from other culture vessels to microscope slides, and thus allows for maintenance of an intact monolayer.

Product Description

Corning BioCoat culture slides, with preapplied extracellular matrix, provide the ability to control *in vitro* cellular environments for cell growth and differentiation under physiologically relevant conditions. They are available with Poly-D-Lysine (PDL), Collagen Type I, PDL/Laminin or Fibronectin.

Falcon culture slides are comprised of a glass slide treated to provide a consistent surface for cell growth and a compartmentalized chamber molded from medical-grade polystyrene. A hydrophobic mask on the slide delineates the location of cells and prevents adhesive from remaining on the culture slides after chamber removal.

A Chamber Removal Key and Safety Removal Fixture are provided to help the user correctly remove the chamber from the glass of both Corning BioCoat and Falcon culture slides.

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ENGLISH

CAUTION

The Removal Key and Safety Removal Fixture **must** be used to remove the chamber from the slide. Please read the removal instructions before using culture slides.

Take precautionary measures during handling; all glass has the potential for breakage. Follow instructions for removal tool and fixture usage.

All samples and collection materials must be handled and disposed of according to the policies and procedures of your facility.

Before the chamber can be removed, culture slides must be exposed for 10 minutes minimum to a 70% to 80% low molecular weight alcohol, such as methanol, ethanol, or isopropanol. Remove chamber immediately. Failure to follow these instructions may result in glass slide breakage.

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ENGLISH

Storage

Store Falcon® culture slides at ambient temperature. Refer to product label for storage temperatures of Corning® BioCoat™ culture slides.

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Chamber Removal Instructions

1. Place the Safety Removal Fixture on a flat surface with the Key Canals and Thumb Tab on top. Fully insert the culture slide, chamber end first, into the Safety Removal Fixture as shown below.

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



Insert
Insérer
Einschieben

ENGLISH

2. Insert the Removal Key prongs into the **key canals**. Once resistance is felt, apply increased pressure until the Removal Key is fully inserted.

3. Lift Chamber from Slide as shown below. Remove Slide and Removal Key.

NOTE: Discard the Removal Key after using three times.

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

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated
Life Sciences

www.corning.com/lifesciences

NORTH AMERICA
t 800.492.1110
t 978.442.2200

ASIA/PACIFIC
Australia/New Zealand
t 61 427286832
Chinese Mainland
t 86 21 3338 4338

India
t 91 124 4604000
Japan
t 81 3-3586 1996
Korea
t 82 2-796-9500
Singapore
t 65 6572-9740
Taiwan
t 886 2-2716-0338

EUROPE
CSEurope@corning.com
France
t 0800 916 882
Germany
t 0800 101 1153
The Netherlands
t 020 655 79 28
United Kingdom
t 0800 376 8660

All Other European Countries
t +31 (0) 206 59 60 51

LATIN AMERICA
grupoLA@corning.com
Brazil
t 55 (11) 3089-7400
Mexico
t (52-81) 8158-8400