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CORNING

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## References

1. Nick JA, et al. Recombinant human activated protein C reduces human endotoxin-induced pulmonary inflammation via inhibition of neutrophil chemotaxis *Blood* 104(13):3878 (2004).
2. Kuijpers TW, et al. Apoptotic neutrophils in the circulation of patients with glycogen storage disease type 1b (GSD1b) *Blood* 101(12):5021 (2003).
3. Kuijpers TW, et al. Neutrophils in Barth syndrome (BTHS) avidly bind annexin-V in the absence of apoptosis *Blood* 103(10):3915 (2004).
4. Slofstra SH, et al. Inhalation of activated protein C inhibits endotoxin-induced pulmonary inflammation in mice independent of neutrophil recruitment *Br J Pharmacol.* 149(6):740 (2006).
5. Stegenga ME, et al. Effect of acute hyperglycaemia and/or hyperinsulinaemia on proinflammatory gene expression, cytokine production and neutrophil function in humans *Diabet Med.* 25(2):157 (2008).
6. Lin Ada, et al. Streptolysin S Inhibits Neutrophil Recruitment during the Early Stages of *Streptococcus pyogenes* Infection and Immunity 77(11):5190 (2009).
7. Johnson Louise A, et al. An inflammation-induced mechanism for leukocyte transmigration across lymphatic vessel endothelium *JEM* 203(12):2763 (2006).

8. Chaudhuri Anathbandhu, STAT1 signaling modulates HIV-1-induced inflammatory responses and leukocyte transmigration across the blood-brain barrier *Blood* 111(4):2062 (2008).
9. Chaudhuri Anathbandhu, 1 HIV-1 activates proinflammatory and interferon-inducible genes in human brain microvascular endothelial cells: putative mechanisms of blood–brain barrier dysfunction *Journal of Cerebral Blood Flow & Metabolism* 28:697 (2008).
10. Lapteva Natalia, et al. Attraction and Activation of Dendritic Cells at the Site of Tumor Elicits Potent Antitumor Immunity *Molecular Therapy* 17(9):1626 (2009).
11. Link Tiffany M, et al. TRPV2 has a pivotal role in macrophage particle binding and phagocytosis *Nature Immunology* 11:232 (2010).
12. Brown MD, et al. Influence of omega-6 PUFA arachidonic acid and bone marrow adipocytes on metastatic spread from prostate cancer *British Journal of Cancer* 102:403 (2010).
13. De Becker Ann, et al. Migration of culture-expanded human mesenchymal stem cells through bone marrow endothelium is regulated by matrix metalloproteinase-2 and tissue inhibitor of metalloproteinase-3 *Haematologica* 92(4):440 (2007).
14. Hassoun Amani T, et al. A rapid screening method for population-specific neuronal motogens, substrates and associated signaling pathways *Journal of Neuroscience Methods* 166:178 (2007).
15. Bernhagen Jürgen, et al. MIF is a noncognate ligand of CXC chemokine receptors in inflammatory and atherogenic cell recruitment *Nature Medicine* 13:587 (2007).
16. Vitari Alberto C, et al. COP1 is a tumour suppressor that causes degradation of ETS transcription factors *Nature* Published online 15 May 2011.
17. Albin A, Noonan DM. The “chemoinvasion” assay, 25 years and still going strong: the use of reconstituted basement membranes to study cell invasion and angiogenesis. *Curr Opin Cell Biol.* 22(5):677-89 (2010).

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**CORNING** | **FALCON**® **AXYGEN**® **PYREX**® **GOSSSELIN**™

For a listing of trademarks, visit [www.corning.com/clstrademarks](http://www.corning.com/clstrademarks).  
All other trademarks are the property of their respective owners.

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[www.corning.com/lifesciences](http://www.corning.com/lifesciences)  
[CLSJP@corning.com](mailto:CLSJP@corning.com)

技術サポート  
Tel:03-3586-1268  
[ScientificSupportJP@corning.com](mailto:ScientificSupportJP@corning.com)

CLS-117-00  
R0-1610-000-T