

# Performance Testing for Axygen® Automation Tip (VTF-384-50UL-R-S)

## Application Note



### Method


The Agilent® Bravo/Velocity11 liquid handling platform was used to assess precision as coefficient of variation (% CV), and accuracy as percent deviation (% D), for Axygen 50 µL tips.

To test the ability of the tip to dispense accurately and precisely at two dispense volumes, 5 µL and 50 µL, a rack of 384 tips aspirated from an Axygen low profile reservoir (Corning Cat. No. RES-SW96-LP) and dispensed into a Corning® 384-well, black, clear bottom microplate (Corning Cat. No. 3711).

For the 5 µL test volume, each tip aspirated 5 µL of Range B solution (Artel Cat. No. MVS-204) and dispensed 5 µL into 50 µL

of diluent solution (Artel Cat. No. MVS-202) in each well. For the 50 µL test volume, each tip aspirated 50 µL of Range A solution (Artel Cat. No. MVS-203) and dispensed 50 µL into 5 µL of diluent solution in each well. To determine the volume of liquid dispensed in each well, absorbance readings for the solutions (diluted Range B solution for 5 µL dispense and Range A solution for 50 µL dispense) were measured using an Artel ELx800NB® plate reader (Artel Cat. No. 1311197). Each study was performed 3 independent times for a total of 1,152 tip dispenses. Evaluation criteria include % D from the set dispense volume and % CV of the measured dispense volume for the 1,152 tip dispenses.

### Results




**Figure 1. Analysis of VTF-384-50UL-R-S tip with aqueous dispense.** The precision (assessed by % CV) and accuracy (assessed by % D) of Axygen VTF-384-50UL-R-S tips dispensing (A) 5 µL and (B) 50 µL volumes using the Agilent Bravo/Velocity11 liquid handling platform were determined using the Artel MVS® system. The % CV and % D were below 2% for both the 5 µL and 50 µL dispenses, n = 1,152.

**Table 1. Aqueous Dispense Results**

Target Volume (µL)	5	50
n	1,152	1,152
% CV	1.78 ± 0.24	1.68 ± 0.29
% D	1.44 ± 0.21	0.15 ± 0.12
Outliers	0	0

### Conclusion

The % CV and % D for the Axygen automation VTF-384-50UL-R-S tips dispensing 5 µL and 50 µL volumes were 5% or below. Therefore, Axygen automation VTF-384-50UL-R-S tips can precisely and accurately dispense volumes as low as 5 µL and as high as 50 µL for aqueous solution using the Agilent Bravo/Velocity11 liquid handling platform.



[www.corning.com/lifesciences/solutions](http://www.corning.com/lifesciences/solutions)

In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs doing what you do every day, across the globe. From collection to analysis, our technical experts understand your challenges and your need for simplified efficient, low- to high-throughput genomics processes.

A combination of global manufacturing expertise, extensive use of in-house automation, an unsurpassed commitment to product innovation and a thorough understanding of your processes enables Corning to offer a beginning-to-end portfolio of high-quality, reliable consumables and reagents for genomics applications.

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences), or contact our Scientific Support Team at [ScientificSupportEMEA@corning.com](mailto:ScientificSupportEMEA@corning.com).

**Corning Incorporated**  
*Life Sciences Europe*

Corning BV  
Fogostraat 12  
1060 LJ Amsterdam  
The Netherlands  
Phone: +31 (0) 20 659 60 51  
Fax: +31 (0) 20 659 76 73  
[CSEurope@corning.com](mailto:CSEurope@corning.com)  
[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

**Support Offices**

**EUROPE**

**France**  
t 0800 916 882  
f 0800 918 636

**Germany**  
t 0800 101 1153  
f 0800 101 2427

**The Netherlands**  
t 31 20 655 79 28  
f 31 20 659 76 73

**United Kingdom**  
t 0800 376 8660  
f 0800 279 1117

**All Other European Countries**  
t 31 (0) 20 659 60 51  
f 31 (0) 20 659 76 73

**CORNING** | **FALCON** | **AXYGEN** | **GOSELIN**

For a listing of trademarks, visit [www.corning.com/clstrademarks](http://www.corning.com/clstrademarks).  
All other trademarks are the property of their respective owners.