

Performance Testing for Axygen® Automation Tip (PK-50-R)

Application Note



Method

The PerkinElmer Janus® liquid handling workstation was used to assess precision as coefficient of variation (% CV), and accuracy as percent deviation (% D), for Axygen 50 µL tips.

To test the ability of the tip to dispense accurately and precisely at two dispense volumes, 5 µL and 50 µL, a rack of 96 tips aspirated from an Axygen low profile reservoir (Corning Cat. No. RES-SW96-LP) and dispensed into a Corning® 96-well, black, clear bottom microplate (Corning Cat. No. 3631).

For the 5 µL test volume, each tip aspirated 5 µL of Range C solution (Artel Cat. No. MVS-205) or DMSO Range C solution

(Artel Cat. No. MVS-217 solution) and dispensed 5 µL into 195 µL of diluent solution (Artel Cat. No. MVS-202) in each well. For the 50 µL test volume, each tip aspirated 50 µL of Range A solution (Artel Cat. No. MVS-203) and dispensed 50 µL into 150 µL of diluent solution in each well. To determine the volume of liquid dispensed in each well, absorbance readings for the solutions (diluted Range C solution for 5 µL dispense and Range A solution for 50 µL dispense) were measured using an Artel ELx800NB® plate reader (Artel Cat. No. 1311197). Each study was performed 3 independent times for a total of 288 tip dispenses. Evaluation criteria include % D from the set dispense volume and % CV of the measured dispense volume for the 288 tip dispenses.

Results




Figure 1. Analysis of PK-50-R tip with aqueous dispense. The precision (assessed by % CV) and accuracy (assessed by % D) of Axygen PK-50-R tips dispensing (A) 5 µL and (B) 50 µL volumes using the PerkinElmer Janus liquid handling workstation were determined using the Artel MVS® system. The % CV and % D were below 1% for both 5 µL and 50 µL dispenses, n = 288.

Table 1. Aqueous Dispense Results

Target Volume (µL)	5	50
n	288	288
% CV	0.80 ± 0.15	0.60 ± 0.14
% D	0.76 ± 0.31	0.87 ± 0.19
Outliers	0	0

Conclusion

The % CV and % D for the Axygen automation PK-50-R tips dispensing 5 µL and 50 µL were 5% or below. Therefore, Axygen automation PK-50-R tips can precisely and accurately dispense volumes as low as 5 µL and as high as 50 µL for aqueous and DMSO solutions using the PerkinElmer Janus liquid handling workstation.





Figure 2. Analysis of PK-50-R tip with DMSO dispense. The precision (assessed by % CV) and accuracy (assessed by % D) of Axygen® PK-50-R tips dispensing 5 µL volumes using the PerkinElmer Janus® liquid handling workstation were determined using the Artel MVS® system. The % CV and % D were below 1% for the 5 µL dispense, n = 288.

Table 2. DMSO Dispense Results

Target Volume (µL)	5
n	288
% CV	0.47 ± 0.03
% D	0.65 ± 0.48
Outliers	0



www.corning.com/lifesciences/solutions

In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs doing what you do every day, across the globe. From collection to analysis, our technical experts understand your challenges and your need for simplified efficient, low- to high-throughput genomics processes.

A combination of global manufacturing expertise, extensive use of in-house automation, an unsurpassed commitment to product innovation and a thorough understanding of your processes enables Corning to offer a beginning-to-end portfolio of high-quality, reliable consumables and reagents for genomics applications.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences, or contact our Scientific Support Team at ScientificSupportEMEA@corning.com.

Corning Incorporated
Life Sciences Europe

Corning BV
Fogostraat 12
1060 LJ Amsterdam
The Netherlands
Phone: +31 (0) 20 659 60 51
Fax: +31 (0) 20 659 76 73
CSEurope@corning.com
www.corning.com/lifesciences

Support Offices

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73



For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.