

# Product Change

---

Date: December 3, 2008  
Subject: Protein Crystallization Microplates

---

**Due to the incompatibility of new solvents used for crystallization, effective immediately Corning Incorporated, Life Sciences will be changing the base resin of the following products from PZero to cyclic olefin (COC). This conversion is expected to be validated by the end of January 2009.**

Corning Cat#	Description
3550	96 Well PZero Protein Crystallization Microplate with 3:1, 1 $\mu$ L Conical Flat Bottom Crystal Cup Wells, Not Treated, Nonsterile
3552	96 Well PZero Protein Crystallization Microplate with 3:1, 2 $\mu$ L Round Bottom Wells, Not Treated, Nonsterile
3553	96 Well PZero Protein Crystallization Microplate with 3:1, 2 $\mu$ L Conical Flat Bottom Wells, Not Treated, Nonsterile