



# Corning Puts Premium Quality Right at Your Fingertips!



## Introducing New and Improved Corning® Centrifuge Tubes with the CentriStar™ Cap

Simply the most innovative, durable and reliable centrifuge tubes in the industry.

### The new CentriStar cap provides...

- ▶ Advanced Ergonomics – wider knurls and roll-over edge design improves grip
- ▶ Leak-Proof Design – revolutionary plug feature that eliminates all seepage
- ▶ Improved Ease-of-Use – quick engagement, easy-on/easy-off, even with one-handed use

### Other advantages include...

- ▶ Highest RCF rating available
  - 15 mL: 12,000 x G
  - 50 mL: 15,500 xG
- ▶ Certified RNase- /DNase-free
- ▶ Greater strength and durability
- ▶ Unique and innovative packaging
  - New universal 15/50 mL rack
  - New bulk, zip lock bags





## Corning® Centrifuge Tubes with CentriStar™ Cap Ordering Information



### 15 mL Conical Centrifuge Tubes, Sterile

Cat. No.	Material	Cap Style	Max RCF (x g)	Inner Pack	Qty/Pk	Qty/Cs
430790	PP	CentriStar Cap	12,000	Foam Rack	50	500
430791	PP	CentriStar Cap	12,000	Zip Lock Bag	25	500



### 50 mL Conical Centrifuge Tubes, Sterile

Cat. No.	Material	Cap Style	Max RCF (x g)	Inner Pack	Qty/Pk	Qty/Cs
430828	PP	CentriStar Cap	15,500	Foam Rack	25	500
430829	PP	CentriStar Cap	15,500	Zip Lock Bag	25	500
4558	PP	CentriStar Cap	15,500	Universal Rack	25	300



### 50 mL Self-Standing Centrifuge Tubes, Sterile

Cat. No.	Material	Cap Style	Max RCF (x g)	Inner Pack	Qty/Pk	Qty/Cs
430921	PP	CentriStar Cap	3,000	Zip Lock Bag	25	500



### Foam Racks

Cat. No.	Description	Qty/Pk	Qty/Cs
431355	Standard Foam Rack for 15 mL Conical Tubes	1	20
4365	Standard Foam Rack for 50 mL Conical Tubes	1	20
4366	Universal Foam Rack for 15 mL and 50 mL Conical Tubes	1	20

## Technical Support

For additional product or technical information, please e-mail us at [CLStechserv@corning.com](mailto:CLStechserv@corning.com), visit our web site [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call 978.635.2200.

## Request a Free Sample

Request a free sample pack online at [www.corning.com/lifesciences/samples](http://www.corning.com/lifesciences/samples) or complete all fields on the form below and fax to **978.635.2484**.

- 15 mL Corning Centrifuge Tube with CentriStar Cap, Polypropylene  
 50 mL Corning Centrifuge Tube with CentriStar Cap, Polypropylene

Name: \_\_\_\_\_

Title: \_\_\_\_\_

Institution: \_\_\_\_\_

Business Address: \_\_\_\_\_

Building/Room/Lab: \_\_\_\_\_

City: \_\_\_\_\_ State/Province: \_\_\_\_\_ Zip/Postal Code: \_\_\_\_\_

Business Telephone: \_\_\_\_\_ E-mail: \_\_\_\_\_

Shipping address if different from above (no P.O. Boxes please): \_\_\_\_\_

City: \_\_\_\_\_ State/Province: \_\_\_\_\_ Zip/Postal Code: \_\_\_\_\_

Primary discipline:  Assay or HTS  BioProcess  Chemistry  Cell/Tissue Culture  Genomics or Microarray

Other \_\_\_\_\_

### Corning Incorporated Life Sciences

45 Nagog Park  
Acton, MA 01720  
U.S.A.

t 800.492.1110  
t 978.635.2200  
f 978.635.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)