

It's not Just a Tube. It's a Lab Legend.

When lab workers are deciding which centrifuge tubes to purchase, they typically prioritize features that impact safety, sample integrity, workflow efficiency, and sustainability.

Here are the most important features, ranked by practical relevance in most lab environments.

FALCON®

A Corning Brand

1.

Leak-proof Seal

No leaks, no screams.

2.

Maximum RCF

Handles up to 16,000 x g* like a roller coaster ride.

*For 50 mL tubes only

3.

Contamination-free

RNase-/DNase-free, nonpyrogenic.

4.

Temperature Resistance

Freezer or heat, I thrive. Storage temperature range is -80°C to 121°C.

5.

Clearly Readable Graduations

You can actually read these markings.

6.

Wide Selection

Available in various capacities.

7.

Environmental Impact

Supply your lab. Sustainably.



Learn more about Corning EcoChoice™ and our centrifuge tubes at ecatalog.corning.com.



www.corning.com/liquidhandling

CORNING

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

The information contained within is accurate as of the date of publication and subject to change without notice. For a listing of trademarks, visit www.corning.com/trademarks. All other trademarks are the property of their respective owners. ©2025 Corning Incorporated. All rights reserved. 10/25 CLS-LH-387 Falcon