Axygen[®] Automation Tips for Hamilton CO-RE[™] Liquid Handling Workstations

96-well Format Tips






1000, 300, 50, and 10 μ L tips

Now there is an alternative source for precision pipet tips compatible with Hamilton automated liquid handling workstations. Corning has a full line of Axygen branded clear and conductive black liquid level sensing tips in 4 sizes: $10~\mu$ L, $50~\mu$ L, $300~\mu$ L, and $1000~\mu$ L to accommodate any scientist's application. Tips are packaged in both hanging tip rack and nested stacked tip racks. As with all Axygen brand robotic tips, they are manufactured to stringent specifications under strict process controls and assembled via automation to ensure consistent performance and quality.

Features and Benefits

- Tips are compatible with Hamilton CO-RE Microlab® Prep™, STAR™, STARlet, STARplus, NIMBUS®, and VANTAGE® workstations
- ▶ The broadest line of tip sizes and styles in the industry
- Unique product choices which are not available from other manufacturers
- Tips engineered to function optimally with no hardware or software adjustments
- ▶ In-house engineering and mold-building/maintenance ensures quality products and a quick customer response
- Manufactured under a Quality Management System which meets the requirements of ISO 9001, which ensures consistent quality and performance
- All tips are validated in-house prior to release ensuring advertised compatibility
- \blacktriangleright 96-well format tips come in 10, 50, 300 and 1,000 μ L sizes (conductive black or clear polypropylene) in hanging or nest stack racks
- ANSI/SLAS footprint rack dimensions for precise fit on the deck
- ▶ With or without aerosol-resistant filters are available
- Nonsterile and sterile versions are available
- ▶ RNase-/DNase-free, human gDNA-free
- Nonpyrogenic






Ordering Information


Products may not be available in all markets.

Axygen® 96-well Format Tips, Hanging Tip Rack

Compatible with Hamilton CO-RE™ Microlab® Prep™, STAR™ Line, Nimbus®, and VANTAGE® Workstations








480

3,840

Conductive Black





Clear

HTF-1000-CBK-HTRS5

Cat. No.	Description	Qty/Pk	Qty/Cs
HT-10-C-HTR-5	10 μL automation tips	480	5,760
HT-10-C-HTR-S-5	10 μL automation tips, sterile	480	5,760
HTF-10-C-HTR-S-5	10 μL filtered automation tips, sterile	480	5,760
HT-50-C-HTR-5	50 μL automation tips	480	5,760
HT-50-C-HTR-S-5	50 μL automation tips, sterile	480	5,760
HTF-50-C-HTR-S-5	50 μL filtered automation tips, sterile	480	5,760


1000 µL liquid level sensing filtered automation tips, sterile


Clear (continued)

Cat. No.	Description	Qty/Pk	Qty/Cs
HT-300-C-HTR-5	300 μL automation tips	480	5,760
HT-300-C-HTR-S-5	300 μL automation tips, sterile	480	5,760
HTF-300-C-HTR-S	300 μ L filtered automation tips, sterile	96	2,304
HTF-300-C-HTR-S-5	300 μL filtered automation tips, sterile	480	5,760
HT-1000-C-HTR-5	1000 μL automation tips	480	3,840
HT-1000-C-HTR-S-5	1000 μL automation tips, sterile	480	3,840
HTF-1000-C-HTR-S	1000 μL filtered automation tips, sterile	96	1,536
HTF-1000-C-HTR-S-5	1000 μL filtered automation tips, sterile	480	3,840

Axygen® 96-well Format Tips, Nested Tip Rack

Compatible with Hamilton CO-RE™ Microlab® Prep™, STAR™ Line, Nimbus®, and VANTAGE® Workstations





Conductive Black

H1
H1
H
Н
Н
Н

Cat. No	Description	Qty/Pk	Qty/Cs
HT-10-CBK-STK	10 μL liquid level sensing automation tips	384	11,520
HT-10-CBK-STK-S	10 μL liquid level sensing automation tips, sterile	384	11,520
HT-50-CBK-STK	50 μL liquid level sensing automation tips	384	11,520
HT-50-CBK-STK-S	50 μL liquid level sensing automation tips, sterile	384	11,520
HT-300-CBK-STK	300 μL liquid level sensing automation tips	384	11,520
HT-300-CBK-STK-S	300 μL liquid level sensing automation tips, sterile	384	11,520

Clear

Cat. No	Description	Qty/Pk	Qty/Cs
HT-10-C-STK	10 μL automation tips	384	11,520
HT-10-C-STK-S	10 μL automation tips, sterile	384	11,520
HT-50-C-STK	50 μL automation tips	384	11,520
HT-50-C-STK-S	50 μL automation tips, sterile	384	11,520
HT-300-C-STK	300 μL automation tips	384	11,520
HT-300-C-STK-S	300 μL automation tips, sterile	384	11,520

Warranty/Disclaimer: Unless otherwise specified, all products are for research use or general laboratory use only.* Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. These products are not intended to mitigate the presence of microorganisms on surfaces or in the environment, where such organisms can be deleterious to humans or the environment. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications. *For a listing of US medical devices, regulatory classifications or specific information on claims, visit www.corning.com/resources.

Corning's products are not specifically designed and tested for diagnostic testing. Many Corning products, though not specific for diagnostic testing, can be used in the workflow and preparation of the test at the customers discretion. Customers may use these products to support their claims. We cannot make any claims or statements that our products are approved for diagnostic testing either directly or indirectly. The customer is responsible for any testing, validation, and/or regulatory submissions that may be required to support the safety and efficacy of their intended application.

CORNING

Corning Incorporated Life Sciences

www.corning.com/lifesciences

NORTH AMERICA t 800.492.1110 t 978.442.2200

ASIA/PACIFIC Australia/New Zealand t 61 427286832 Chinese Mainland t 86 21 3338 4338 India t 91 124 4604000 Japan t 81 3-3586 1996

Korea t 82 2-796-9500 Singapore t 65 6572-9740 Taiwan

t 886 2-2716-0338

CSEurope@corning.com
France
t 0800 916 882
Germany
t 0800 101 1153
The Netherlands

EUROPE

t 020 655 79 28 **United Kingdom** t 0800 376 8660 All Other European Countries t+31 (0) 206 59 60 51

LATIN AMERICA grupoLA@corning.com

Brazil t 55 (11) 3089-7400 Mexico

t (52-81) 8158-8400

Hamilton®, Microlab®, NIMBUS®, and VANTAGE® are registered trademarks of Hamilton Company. CO-RE™, Prep™, STAR™, STARIet™, and STARplus™ are trademarks of Hamilton Company. The compatibility of Axygen tip products with Hamilton Company products has been independently determined by Corning; however, there should be no implication of sponsorship or endorsement by Hamilton Company of Axygen tip products.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners. © 2012-2024 Corning Incorporated. All rights reserved. 11/24 CLS-A-AUTO-013 REV10