

Corning® glutagro™ Supplement

CORNING



The Best Research Requires The Healthiest Cells.

Corning glutagro supplement combines the quality you expect from Corning with the reagent stability you need. Corning glutagro supplement is the perfect replacement for L-glutamine.

Cells grown in traditional media supplemented with L-glutamine are susceptible to cytotoxicity. In media supplemented with Corning glutagro, cells receive this needed supplement in the healthiest environment possible.


This modified form of L-glutamine does not degrade into ammonia during storage and incubation, minimizing cytotoxicity, and improving culture quality.

Ordering Information

Cat. No.	Description	Unit Size	Qty/Pk
25-015-CI	Corning glutagro supplement	100 mL	1
10-101-CV	DMEM with Corning glutagro	500 mL	6
10-102-CV	DMEM with Corning glutagro, without sodium pyruvate	500 mL	6
10-103-CV	DMEM/F-12 with Corning glutagro	500 mL	6
10-104-CV	RPMI 1640 with Corning glutagro	500 mL	6

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.



www.corning.com/lifesciences/solutions

At Corning, cells are in our culture. In our continuous efforts to improve efficiencies and develop new tools and technologies for life science researchers, we have scientists working in Corning R&D labs across the globe, doing what you do every day. From seeding starter cultures to expanding cells for assays, our technical experts understand your challenges and your increased need for more reliable cells and cellular material.

It is this expertise, plus a 160-year history of Corning innovation and manufacturing excellence, that puts us in a unique position to offer a beginning-to-end portfolio of high-quality, reliable cell culture consumables.

For additional product or technical information, visit www.corning.com/lifesciences or call 1.800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

Mediatech, Inc.
A Corning Subsidiary
9345 Discovery Boulevard
Manassas, VA 20109
t 800.235.5476
t 703.471.5955
f 703.467.9851
www.corning.com/lifesciences/media

CORNING | **FALCON** | **AXYGEN** | **GOSSELIN** | **PYREX**

For a listing of trademarks, visit www.corning.com/clstrademarks.
All other trademarks are the property of their respective owners.