

December 12, 2024

To Whom It May Concern,

As an accommodation to our customers, Corning Life Sciences has performed testing and verification to establish a *minimum* shelf life for certain Axygen®, Corning®, Costar®, Gosselin® and Falcon® products. These products have not been tested past this established minimum shelf life. The customer is advised to evaluate the product for conformity with specifications and suitability for the use intended if there is any intention to use the product past the established minimum shelf life.

Shelf life is based on the date of manufacture. The manufacture date can be found on the packaging label and on the certificate of conformance. If the product is also labeled, the product label has a date code on it. The expiration date is calculated based on the manufacture date and shelf life and is printed on the certificate of compliance, secondary packaging and in some instances primary packaging.

The established shelf life is only applicable to contents of unopened and undamaged packaging. This cautionary statement is also printed on the sterile barrier (primary packaging) of most products. We recommend that products be stored in ambient conditions (unless otherwise indicated), protected from humidity, exposure to UV light, ozone, exhaust fumes and paint fumes.

After expiration of the *minimum* shelf life, or in situations where the product may not have been handled or stored properly, the customer is advised to evaluate the product for conformity with specifications and suitability for the use intended.

If you have any questions or concerns, please contact us at
ScientificSupport@corning.com

Thank you for your business with Corning Life Sciences.



Neil Van Os
Senior Life Sciences Quality Assurance Manager,
Commercial & Compliance