
# Corning<sup>®</sup> Custom Media and Serum Solutions

Media, sera, and reagents to meet your unique needs from research to commercial production

CORNING



# Corning® Custom Media and Serum Solutions

Leverage Corning's expertise and broad product portfolio to support your product pipeline from discovery to commercial production. From cell culture and drug discovery to vaccine and antibody production, our custom manufacturing services can produce tailored media, sera, reagents, and sterile solutions to meet your unique research and production needs.

Media products are manufactured under the current ISO 13485 standard and FDA Quality System Regulation 21 CFR 820, also known as cGMP (current good manufacturing practices). We offer industry-leading media products and process quality that meets cGMP requirements for manufacturing and sterility.

Our serum products are collected, processed, and tested in accordance with cGMP standards and EP (European Pharmacopeia) and USP (United States Pharmacopeia) methods. Corning facilities are ISIA traceability certified for serum. Some serum products have Certificates of Suitability from EDOM.

Corning custom media and sterile solutions are produced in separate designated animal and non-animal manufacturing suites and high quality water is used in all production. Corning custom media and sterile solutions can be supplied in liquid or powder formats. Our expert manufacturing and reliable supply of serum products perfectly complement our classical cell culture media portfolio.

# Large-scale cGMP Batch Capabilities

Packaging/Container Sizes	Volume
Liquid Media	
Overall batch sizes	100L to 8,000L
Animal	10L to 200L 1,000L to 6,000L
Non-animal	100L to 200L 400L to 2,000L 1,600L to 8,000L
Powder Media	
Overall batch sizes	50 to 1,000 kg

#### **Custom Small-scale Pilot Batches**

We can rapidly provide small-scale pilot batches ideal for development and testing purposes prior to cGMP scale-up. A summary of quality control tests performed on liquid pilot batches is available for endotoxin, appearance, pH, osmolality, and conductivity.

Packaging/Container Sizes	Volume
Standard Bags	
1L, 5L, 10L, 20L bags	1L to 100L
Liquid Media	
Standard bottles: 50 mL, 100 mL, 500 mL, 1L bottles	1L to 100L





Corning can provide you with a full range of products to meet your custom media formulation, packaging, and regulatory requirements and support you with experience and expertise along the way.

# Our standard and custom product offerings include:

## **Media and Reagents**

- Classical media and serum-free formulations
- Sterility challenge media (tryptic soy broth) for use in biopharmaceutical production validation
- Supplements
- Insect cell culture formulations
- Downstream purification buffers
- Cell and tissue processing solutions
- Diluents and controls
- Customer specific formulations

#### **High Quality Water**

- Tested to specifications that may include USP and/or EP and/or JP monograph requirements
- Cell culture grade
- Molecular grade: RNase-/DNase-free and protease-free
- ▶ Reagent grade: RNase-/DNase-free and protease-free

# Our serum product offerings include:

## FBS (Fetal Bovine Serum)

- Corning's premium FBS is only collected from government approved abattoirs in the United States, Australia, and New Zealand.
- Corning's regular FBS is derived from government registered facilities in North and Central America.
- Traceability certified according to the ISIA Traceability Program.
- Certificate of Suitability (CEP) issued for United States origin.
- Comprehensive testing meeting USP and EP requirements and flexibility to complete custom needs.
- ▶ Tank heat inactivation technology.
- ▶ Validated gamma irradiation protocols in place for two dose ranges of 25 to 40 kGy and 30 to 45 kGy.
- Extensive portfolio of treatments available including dialysis, charcoal dextrane, and IgG removal.

#### **Bovine Calf Serum**

- Iron fortification capability upon request.
- ▶ Tank heat inactivation technology.
- Validated gamma irradiation protocols in place for two dose ranges of 25 to 40 kGy and 30 to 45 kGy.

In addition to FBS and Bovine Calf Serum, Donor Horse, Sheep, and Rabbit Serum are available upon request.

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# **CORNING**

**Corning Incorporated** Life Sciences

836 North St. Building 300, Suite 3401 Tewksbury, MA 01876 t 800.492.1110 t 978.442.2200 f 978.442.2476 www.corning.com/lifesciences

ASIA/PACIFIC Australia/New Zealand t 61 427286832 Chinese Mainland t 86 21 3338 4338 f 86 21 3338 4300 India t 91 124 4604000 f 91 124 4604099

Japan t 81 3-3586 1996 f 81 3-3586 1291 Korea t 82 2-796-9500 f 82 2-796-9300 Singapore t 65 6572-9740 f 65 6735-2913 Taiwan t 886 2-2716-0338 f 886 2-2516-7500 EUROPE CSEurope@corning.com France t 0800 916 882 f 0800 918 636 Germany t 0800 101 1153 f 0800 101 2427 The Netherlands t 020 655 79 28

t 0800 376 8660 f 0800 279 1117

f 020 659 76 73 **United Kingdom**  All Other European Countries t +31 (0) 206 59 60 51 f+31 (0) 206 59 76 73

LATIN AMERICA grupoLA@corning.com Brazil t 55 (11) 3089-7400 Mexico t (52-81) 8158-8400