

CORNING

Meet Flaska:
The Flask That
Gets You.

I'm not just a flask.
I'm your lab partner!



Another Day, Another Spill

Why is working
with flasks a three-
handed job?



Sounds like
someone needs
an upgrade.

Flaska, the U-shaped Hero



I'm Flaska – here to
keep you and your
cells happy.



- Rounded shoulders**
Easier rinsing and
cells scraping
- Shorter height**
Better incubator
stacking
- Wider mouth**
No more scrapers or
pipet gymnastics

Less Struggle. More Science.



- Large clear viewing surface
- Canted neck for smooth pouring
- One-handed media handling
- Clear graduations
- Easy cap access for removal and tightening

I told you – better
access, better results,
better day!

Flaska's Secret Formula

- Made from virgin USP Class VI polystyrene
- Nonpyrogenic and BSE/TSE-certified
- Sterility Assurance Level (SAL) of 10⁻⁶
- Manufactured under cGMP standards
- Available in TC-treated, non-treated, and Corning CellBIND® surfaces



Flaska's Green Advantage



Same surface area.
23% less plastic.*
No compromise.

- Sustainable switch = less waste,
more room for your science.
- Stable stacking = fewer flask-towers
in the incubator.

*U-shaped design with 75 cm² uses 23% less plastic than our standard T-75 flasks, and U-shaped design with 150 cm² uses 30% less plastic than our standard T-150 flasks – while keeping the same cell culture surface area.

Join the U-Revolution. Level up with Flaska. Your cells will thank you.



www.corning.com/lifesciences

CORNING