

オートメーションチップ互換性チャート

ワークステーション ブランド	モデル	EV-50 EVF-50	EV-100 EVF-100	EV-200 EVF-180	EVF-205-WB	FDSS-1536	FX-10	FX-20 FXF-20	FX-50 FXF-50	FX-250 FXF-180 FXF-200	FX-255-WB FXF-165-WB	FX-384 FXF-384	FX-384-XL FXF-384-XL	FX-384-30FP FXF-384-30FP	FX-384-45FP-R	FX-1536-30FP	HT-10 HTF-10	HT-50 HTF-50	HT-300 HTF-300	HT-1000 HTF-1000	
Aurora Biomed	VERSA™																				
Agilent/Velocity11	Bravo VPrep®																				
Beckman Coulter	Biomek® 1,000																				
	Biomek 2000																				
	Biomek 3000						●	●	●	●	●										
	Biomek 4000						●	●	●	●	●										
	Biomek FX-96/NX-96						●	●	●	●	●										
	Biomek FX-384/NX-384											●	●	●	●	●					
	Multimek™ AP-96						●	●	●	●	●										
	Multimek AP-384											●	●	●	●	●					
	Multimek CCS-head Multimek Type II																				
BioTek	Precision™	●	●	●	●																
	Precision XS	●	●	●	●																
Corbett	CAS-1200																				
	CAS-4200, CAS-4800																				
Dynamic Devices	Oasis LM	●	●	●	●							●	●	●	●						
Hamamatsu	7000					●															
Molecular Devices	FLIPR® Liberty 384																				
	FLIPR Tetra																				
Hamilton	Nimbus®											●	●	●	●		●	●	●	●	
	Microlab® STAR™											●	●	●	●		●	●	●	●	
	Microlab STARlet											●	●	●	●		●	●	●	●	
	Microlab STARplus											●	●	●	●		●	●	●	●	
PerkinElmer	Evolution™ P3																				
	Janus®																				
	Multiprobe II ET (hanging wafer)																				
	Multiprobe II HT																				
	Platetrak™																				
	Zymark® RapidPlate®	●	●	●	●																
	Sciclone® Zephyr®	●	●	●	●																
ProGroup/Wellpro	3000-96											●	●	●	●						
	3000-384											●	●	●	●						
Tecan	Aquarius	●	●	●	●																
	Freedom EVO®	●	●	●	●																
	Genesis Freedom®																				
	Miniprep with LiHa																				
Te-Mo	●	●	●	●																	
Qiagen	BioRobot®																				

