


## Starting A New Lab?

# Get 25% Free Goods and Technical Training



Get your lab started with Corning®, Costar® and PYREX® products and we will give you free goods equal to 25% of your total purchase plus a CD-ROM containing cell culture training courses.\*

### Cell Culture Training Courses on CD

Corning has developed these cell culture training courses by combining over 30 of the most popular Corning articles, protocols, recorded seminars and product selection guides. These are designed to help researchers, technicians and students learn more about the skills and techniques required to successfully use cell culture in the laboratory. A series of optional study questions is provided to help determine if students or employees have mastered the key concepts covered in these articles.

### Example of 25% Free Goods Offer

Spend \$1,000 on Corning and Costar products and get \$250 in free Corning and Costar products.

### To Claim Your Free Product:

- ▶ **By E-mail:** Include the below information along with an electronic copy of your proof of purchase (invoice, order confirmation, etc.) in an e-mail and send to [Lspromo@corning.com](mailto:Lspromo@corning.com). Include "New Lab-GWEB" in the subject line.
- ▶ **By Fax:** Complete the below information and fax it with a copy of your proof of purchase to **978.442.2484**.

Name: \_\_\_\_\_ Title: \_\_\_\_\_

Institution/Company: \_\_\_\_\_

Business address (No P.O. boxes please): \_\_\_\_\_ Bldg./Room/Lab: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Business telephone: \_\_\_\_\_ E-mail: \_\_\_\_\_

Shipping address (if different than business address): \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

**Primary Activity:**  Assay or HTS  Bioprocess  Cell/Tissue Culture  Chemistry or Biochemistry  Genomics or Microarrays  
 Other \_\_\_\_\_

Please send me the Cell Culture Training CD.

Free Goods: Select from checklist on back or from any qualified product\* in Corning's catalog.

Qty.: \_\_\_\_\_ Cat. No.: \_\_\_\_\_ Qty.: \_\_\_\_\_ Cat. No.: \_\_\_\_\_

Qty.: \_\_\_\_\_ Cat. No.: \_\_\_\_\_ Qty.: \_\_\_\_\_ Cat. No.: \_\_\_\_\_

Qty.: \_\_\_\_\_ Cat. No.: \_\_\_\_\_ Qty.: \_\_\_\_\_ Cat. No.: \_\_\_\_\_

Please indicate which catalog you are using to reference the catalog numbers above:

Corning Catalog  Dealer Catalog (indicate dealer): \_\_\_\_\_

P.O. #(s): \_\_\_\_\_

For additional product information or technical documents offered by Corning, please visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

\*Offer is good for all new labs for their initial stocking order and up to 2 subsequent orders within 3 months of the original order. Only one CD will be issued per lab. Minimum order value of \$2,000 applies to CD only. The free goods provided will be your choice of any Corning and Costar plastic disposable or glass labware. Qualifying Corning and Costar products include all disposable lab products, tissue culture plastics, PYREX® glassware, hot plates and pipettors.

CORNING

Choose your free product from these or any qualified product\* in our catalog.

Cat. No.	Description	Qty/Cs
<b>Corning® Tissue Culture Flasks</b>		
430639	TC Flask 25 cm <sup>2</sup> Canted Neck, Anti-Tip, Vent, Sterile	200
430641	TC Flask 75 cm <sup>2</sup> Canted Neck, Anti-Tip, Vent, Sterile	100
430823	TC Flask 150 cm <sup>2</sup> Canted Neck, Anti-Tip, Plug Seal, Sterile	50
<b>Costar® Multiple Well Plates</b>		
3516	Multiple Well Cluster Plate, 6 Well, TC-Treated, Sterile	50
3513	Multiple Well Cluster Plate, 12 Well, TC-Treated, Sterile	50
3524	Multiple Well Cluster Plate, 24 Well, TC-Treated, Sterile	100
3548	Multiple Well Cluster Plate, 48 Well, TC-Treated, Sterile	100
3596	Multiple Well Cluster Plate, 96 Well, TC-Treated, Sterile, Individual	50
<b>Corning Cell Culture Dishes</b>		
430166	Cell Culture Dish, 60 mm Style, TC Treated, Sterile	20
430167	Cell Culture Dish, 100 mm Style, TC Treated, Sterile	20
<b>Corning Stripette® Serological Pipets</b>		
4487	Stripette Pipet, Individual 5 mL, Paper/Plastic Wrap	200
4488	Stripette Pipet, Individual 10 mL, Paper/Plastic Wrap	200
4489	Stripette Pipet, Individual 25 mL, Paper/Plastic Wrap	200
4490	Stripette Pipet, Individual 50 mL, Paper/Plastic Wrap	100
<b>Corning 15 mL Centrifuge Tubes</b>		
430791	Centrifuge Tube 15 mL, PP, Flat Top, Sterile, Bulk, 12,000 RCF	500
<b>Corning 50 mL Centrifuge Tubes</b>		
4558	Centrifuge Tube, 50 mL, PP, CentriStar™ Cap, Universal Rack, 15,500 RCF	300
430829	Centrifuge Tube, 50 mL, PP, CentriStar Cap, Bulk, 15,500 RCF	500
<b>Corning Bottle Filter Systems</b>		
431096	Filter System, 250 mL, 0.22 μm, PES	12
431097	Filter System, 500 mL, 0.22 μm, PES	12
431098	Filter System, 1000 mL, 0.22 μm, PES	12
<b>Corning Cryogenic Vials</b>		
430487	Cryogenic Vial, Int-Thread, 1.2 mL, Round, Self-Standing	500
430659	Cryogenic Vial, Ext-Thread, 2.0 mL, Round, Self-Standing	500
430658	Cryogenic Vial, Ext-Thread, 1.2 mL, Round, Self-Standing	500
430488	Cryogenic Vial, Int-Thread, 2.0 mL, Round, Self-Standing	500
<b>Corning Microcentrifuge Tubes</b>		
3620	Microcentrifuge Tube, 1.7 mL, Snap Cap, Clear	500
3622	Microcentrifuge Tube, 1.7 mL, Snap Cap, Rainbow	500
<b>Corning Hinged Racked Pipet Tips</b>		
4711	Pipet Tip, 1-200 μL, 96/Hinged Rack, Yellow, Sterile	960
4714	Pipet Tip, 100-1000 μL, 100 Tip, Hinged Rack, Sterile	1000
<b>Corning Non-Hinged Racked Pipet Tips</b>		
4894	Pipet Tip, 0.1-10 μL, Non-hinged Rack, Natural, Sterile	960
4864	Pipet Tip, 1-200 μL, Non-hinged Rack, Natural, Sterile	960
9032	Pipet Tip, 100-1000 μL, 100 Tip, Non-hinged Rack, Sterile	1000
<b>Corning Filter Pipet Tips</b>		
4807	IsoTip, 0.1-10 μL, G Tip, 96/Rack, Sterile, Gilson/Rainin Style Fit	960
4821	IsoTip, 1-30 μL 96/Rack, Sterile, Universal Fit	960
4823	IsoTip, 1-200 μL 96/Rack, Sterile, Universal Fit	960
4809	IsoTip, 100-1000 μL 100/Rack, Sterile, Universal Fit	1000
<b>PYREX® Round Media/Storage Bottles</b>		
1395-100	100 mL Round Storage Bottle, Screw Cap	10
1395-250	250 mL Round Storage Bottle, Screw Cap	10
1395-500	500 mL Round Storage Bottle, Screw Cap	10
1395-1L	1000 mL Round Storage Bottle, Screw Cap	10
1395-2L	2000 mL Round Storage Bottle, Screw Cap	10
<b>PYREX Beaker, Low Form Double Scale Graduated</b>		
1000-Pack	Beaker, Assortment Pack, 1 each, 50 mL, 100 mL, 250 mL, 600 mL, 1L	1 set
1000-100	100 mL Beaker	48
1000-250	250 mL Beaker	48
1000-600	600 mL Beaker	36
1000-1L	1L Beaker	24
<b>PYREX Graduated Cylinders</b>		
3024-Pack	Cylinder, Assortment Pack, 1 each, 10 mL, 25 mL, 50 mL, 100 mL, 250 mL	1 set
3024-25	25 mL Cylinder	18
3024-50	50 mL Cylinder	18
3024-100	100 mL Cylinder	12
3024-250	250 mL Cylinder	12
<b>PYREX Erlenmeyer Flasks, Narrow Mouth, Heavy Duty Rim, Graduated</b>		
4980-Pack	Erlenmeyer Flask Assortment Pack – 1each: 50 mL, 125 mL, 250 mL, 500 mL, 1L	1 set
4980-125	125 mL Erlenmeyer Flask	48
4980-250	250 mL Erlenmeyer Flask	8
4980-500	500 mL Erlenmeyer Flask	36
4980-1L	1000 mL Erlenmeyer Flask	24
<b>Corning 96 Well Microplate</b>		
3603	Assay Plate, 96 Well, Black with Clear Flat Bottom, Tissue Culture Treated, Sterile, with Lid	1 pack

Visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) for current prices.

Corning, Costar, PYREX, and Stripettes are registered trademarks of Corning Incorporated, Corning, NY. CentriStar, HYPERFlask, and RoboFlask are trademarks of Corning Incorporated, Corning, NY.

Corning Incorporated, One Riverfront Plaza, Corning, NY 14831-0001

**CORNING**

Corning Incorporated  
Life Sciences

Tower 2, 4th Floor  
900 Chelmsford St.  
Lowell, MA 01851  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)